

24-Hour Emergency Number
(ChemTel) 888.255.3924 in U.S.

cleaning@alconox.com

alconox.com

- distributors
- technical information
- free samples
- new developments
- white papers
- prices



TERGAZYME® Enzyme-Active Powdered Detergent

- Concentrated to save you money
- Biodegradable and readily disposable
- Replaces corrosive acids and hazardous solvents
- Protease enzyme removes proteinaceous soils, tissue, blood and body fluids
- Free rinsing to provide reliable results and no interfering residues
- Use to pass cleaning validation tests for lab accreditation and plant inspection approval

Used to clean: Hospital instruments, dairy equipment, laboratory ware, reverse osmosis and ultrafiltration membranes and units, sampling apparatus, pharmaceutical apparatus, cosmetics manufacturing equipment, tubing, pipes, optical parts, process equipment, industrial parts, desalination plants, tanks and reactors. Authorized by USDA for use in federally inspected meat and poultry plants. Passes inhibitory residue test for water analysis. FDA certified.

Used to remove: Soil, grit, grime, blood, tissue, grease, fats, oils, proteinaceous soils, dairy proteins, particulates, solvents and bioreactor residue.

Surfaces cleaned: Corrosion inhibited formulation recommended for glass, metal, stainless steel, porcelain, ceramic, plastic, rubber and fiberglass. Can be used on soft metals such as copper, aluminum, zinc and magnesium if rinsed promptly. Corrosion testing may be advisable.

Cleaning method: Soak, brush, sponge, cloth, ultrasonic, flow through clean-in-place. Will foam—not for spray or machine use.

Directions: Make a fresh 1% solution (2 1/2 tbsp per gal, 1 1/4 oz per gal, or 10 g per L) in cold or warm water. If available, use warm water below 130°F (55°C). Clean by soak, circulate, wipe, or ultrasonic method. Follow manufacturers directions for filter membrane cleaning. Not for spray machines: will foam. RINSE THOROUGHLY—preferably with running water. For critical cleaning, do final or all rinsing in distilled,

deionized, or purified water. For food contact surfaces, rinse with potable water. Used on a wide range of glass, ceramic, plastic, and metal surfaces. Corrosion testing may be advisable. Material compatibility guidelines available on request.

Convenient Sizes	Alconox Cat. #
Case 9 x 4 lb (9 x 1.8 kg) boxes	1304
25 lb (11 kg) carton	1325
50 lb (22 kg) carton	1350
100 lb (45 kg) drum	1301
300 lb (136 kg) drum	1303
250 x 1/2 oz box	1312
1 lb makes 13 gal (1 kg makes 100L) cleaning solution	



TERGAZYME is available from leading laboratory, hospital, clinical and industrial suppliers. To find a distributor for Alconox, Inc. detergents, use the Locate Dealers search tool in the Support section of alconox.com. To request FREE samples, use the Get Samples form at alconox.com. Or call Alconox, Inc. at +914.948.4040, today.

TAZTB.1

24-Hour Emergency Number
(ChemTel) 888.255.3924 in U.S.
cleaning@alconox.com

alconox.com

- distributors
- technical information
- free samples
- new developments
- white papers
- prices

Physical Data	Typical Value
pH of 1% Solution	9.5
Surface Tension 1% Sol'n (dyne/cm)	32
Conductivity 1% Sol'n (mS/cm at 21 degrees C)	6.3
Form	Powder
Color	White and cream colored flakes and brown specks
Solubility in Water	To 10% (w/w) at ambient temperature
Hard Water Effectiveness	Highly effective
Biodegradability	Yes, meets EC648/2004
Foam Tendency	High foaming
Shelf Life	24 months from end of month of manufacture
Fragrance Content	0%
Phosphate Content (as phosphorus)	7.5%
Flash Point (degrees F)	None
Organic Carbon (% calculated w/w)	11%

Chemical Description

TERGAZYME consists primarily of a homogeneous blend of sodium linear alkylaryl sulfonate, phosphates, carbonates, and protease enzyme. TERGAZYME is anionic in nature. The protease enzyme in TERGAZYME is Bacillus licheniformis subtilisin Carlsberg which may be deactivated by 300 ppm hypochlorite at 85 degrees F in seconds; 3.5 ppm hypochlorite at 100 degrees F for 2 min; exposure to pH below 4 for 30 min at 140 degrees F; or by heating to 175 degrees F for 10 min.

Cleaning Validation Methods

To quantify trace levels of Alconox for cleaning validation purposes, use conductivity, UV, TOC, HPLC, or phosphorus detection. For details see the Cleaning Validation page in Support at alconox.com. For lab accreditation tests to show the detergent is rinsed away, see the Accreditation Residue Detection in the Support section of alconox.com.

Health Safety Information

OSHA Hazardous Ingredients	None
RCRA Hazard Class	Nonhazardous
EPA Priority Pollutants	None
DOT Hazard Class	Nonhazardous
Flammability	Nonflammable
Latex Content	None in product or packaging.
Irritation	Cat. 2 skin, Cat. 2A eye irritant
Oral Toxicity	LD ₅₀ > 5000 mg/kg oral rat.
Inhalation Toxicity	Solution: No, Powder dust: irritant
VOC Content	0%
GMO Content	None
BSE/TSE Source	None
2007/68/EC Annex IIIa Allergens	None
REACH Art. 57 (SVHC)	None
RoHS	No RoHS hazardous ingredients.
Carcinogenicity	NTP = No IARC = No OSHA = No
Ingredients on TSCA Inventory	Yes
EC 907/2006	Contains no enzymes, optical brighteners, perfumes, allergenic fragrances, or preservative agents.

Precautions

No special precautions other than good industrial hygiene and safety practices employed with any industrial chemical (see Directions). A Safety Data Sheet is available at alconox.com.

Contact Alconox, Inc. for purchase specifications. Information presented is typical. Not to be taken as specifications. Typical data is not a specification.

While the information in this report should not be considered to be a product warranty, we urge you to investigate, test and verify the suitability of Alconox, Inc. detergents for your specific application. We, of course, can not give permission to use, or recommend the use of, our detergents where they infringe patents. No representation or warranty is made as to the safety of products or materials mentioned under the Federal Food Additives Amendment of 1958.

TAZTB.1