**Thermo Scientific  
Laboratory  
Temperature  
Control Products**

Manual Part Number U01361  
Rev. 03/03/2021

**ADVANCED Series Heated  
Immersion Circulators**

AC150      AC200

**ARCTIC Series  
Refrigerated/Heated  
Bath Circulators**

A10            A25            A45HC

A10B          A25B

**SAHARA Series Heated  
Bath Circulators**

S7            S13            S21

S45          S49

**Installation**

**Operation**

**Basic Maintenance**

Visit our Web site at:

<http://www.thermofisher.com/tc>  
Product Service Information, Applications  
Notes, SDS Forms, e-mail.

**Thermo Scientific  
Laboratory  
Temperature  
Control Products**

Manual Part Number U01361  
Rev. 03/03/2021

**ADVANCED Series Heated  
Immersion Circulators**

AC150      AC200

**ARCTIC Series  
Refrigerated/Heated  
Bath Circulators**

A10            A25            A45HC  
A10B          A25B

**SAHARA Series Heated  
Bath Circulators**

S7            S13            S21  
S45          S49

**Installation  
Operation  
Basic Maintenance**

Visit our Web site at:

<http://www.thermofisher.com/tc>  
Product Service Information, Applications  
Notes, SDS Forms, e-mail.

## **Thermo Fisher Scientific**

25 Nimble Hill Road  
Newington, NH 03801  
Tel : (800) 258-0830 or  
(603) 436-9444  
Fax : (603) 436-8411  
[www.thermoscientific.com/tc](http://www.thermoscientific.com/tc)

## **Sales, Service, and Customer Support**

25 Nimble Hill Road  
Newington, NH 03801  
Tel: (800) 258-0830  
Sales: 8:00 am to 5:00 pm  
Service and Support: 8:00 am to 6:00 pm  
Monday through Friday (Eastern Time)  
Fax: (603) 436-8411  
[service.tc.us@thermofisher.com](mailto:service.tc.us@thermofisher.com)

Dieselstrasse 4  
D-76227 Karlsruhe, Germany  
Tel : +49 (0) 721 4094 444  
Fax : +49 (0) 721 4094 300  
[info.tc.de@thermofisher.com](mailto:info.tc.de@thermofisher.com)

Building 6, No. 27  
Xin Jinqiao Rd., Shanghai 201206  
Tel : +86(21) 68654588  
Fax : +86(21) 64457830  
[info.china@thermofisher.com](mailto:info.china@thermofisher.com)

### **Statement of Copyright**

Copyright © 2021 Thermo Fisher Scientific. All rights reserved.

This manual is copyrighted by Thermo Fisher Scientific.

Users are forbidden to reproduce, republish, redistribute, or resell any materials from this manual in either machine-readable form or any other form.

# Table of Contents

<b>Preface</b>	<b>i</b>
Compliance .....	i
WEEE .....	i
After-Sale Support .....	ii
Feedback.....	ii
Unpacking .....	ii
<b>Section 1 Safety</b>	<b>1-1</b>
Safety Warnings .....	1-1
<b>Section 2 General Information</b>	<b>2-1</b>
Description and Intended Use.....	2-1
Bath/Circulator Specifications.....	2-1
Wetted Materials.....	2-8
<b>Section 3 Installation</b>	<b>3-1</b>
Ambient Conditions .....	3-1
Ventilation .....	3-1
Bath Circulator .....	3-2
Electrical Requirements .....	3-3
Grounding Strap and Nut .....	3-3
External Circulation .....	3-7
Remote Temperature Sensor.....	3-8
Multifunction Port (AC 200 only) .....	3-8
USB Port (AC 200 only) .....	3-8
Tubing Requirements .....	3-8
Tubing.....	3-9
Fluids .....	3-10
Additional Fluid Precautions .....	3-12
Filling Requirements .....	3-15
Draining.....	3-15
<b>Section 4 Operation</b>	<b>4-1</b>
ADVANCED Heater Immersion Circulator .....	4-1
Setup .....	4-2
Start Up .....	4-2
Status Display .....	4-3
Stand By Mode.....	4-3
Stopping the Circulator.....	4-4
Power Down.....	4-4
Shut Down.....	4-4
Restarting .....	4-4

	Changing the Setpoint.....	4-5
	Menu Displays.....	4-6
	Menu.....	4-6
	Menu Tree.....	4-7
	Settings - Application. Settings.....	4-8
	Alarms.....	4-9
	Fluid Type.....	4-10
	Pump Adjustment.....	4-10
	Set Clock.....	4-11
	On/Off Timer.....	4-11
	Remote Sensor.....	4-11
	Auto Restart.....	4-11
	Energy Savings.....	4-11
	Edit - Ramp Program.....	4-12
	Running a Ramp Program.....	4-14
	Settings - Display Options.....	4-16
	System - Messages.....	4-17
	System - Configuration.....	4-17
	System - Password Reset.....	4-18
	System - Accessory.....	4-21
	High Temperature Cutout.....	4-22
	Decommissioning/Disposal.....	4-23
	Storage.....	4-23
<b>Section 5</b>	<b>Preventive Maintenance.....</b>	<b>5-1</b>
	Cleaning.....	5-1
	Condenser Fins.....	5-1
	Tubing.....	5-1
	Grounding Strap and Nut.....	5-2
	Testing the Safety Features.....	5-2
<b>Section 6</b>	<b>Troubleshooting.....</b>	<b>6-1</b>
	Error Displays.....	6-1
	Checklist.....	6-5
<b>Appendix</b>	<b>AC Serial Communications.....</b>	<b>A-1</b>
<b>Warranty</b>		

# Preface

## Compliance

### US/Canada

UL Listed to US and Canadian Standards.

**Note** The cULus Mark on the nameplate on the bath applies to the control head and bath covered by the BOM number on the nameplate, see below. The cULus Mark does not apply to the bath or control head when sold separately, or if the user changes one or the other. ▲

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### European Union - CE

The Declaration of Conformity is located in the back of this manual.

## WEEE

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2012/19/EU. It is marked with 'wheelie bin' symbol:

Thermo Fisher Scientific has contracted with one or more recycling/disposal companies in each EU Member State, dispose of or recycle this product through them. Further information on Thermo Fisher Scientific's compliance with these Directives is available at:

[www.thermoscientific.com/WEEERoHS](http://www.thermoscientific.com/WEEERoHS)

## After-sale Support

Thermo Fisher Scientific is committed to customer service both during and after the sale. If you have questions concerning the bath operation, or questions concerning spare parts or Service Contracts, call our Sales, Service and Customer Support phone number, see this manual's inside cover for contact information.

There can be up to three nameplates located on the rear of the bath. Before calling, please obtain the serial number printed on the complete system nameplate located on the upper rear of the bath.

Sample Nameplate

## Feedback

We appreciate any feedback you can give us on this manual. Please e-mail us at [tcmanuals@thermofisher.com](mailto:tcmanuals@thermofisher.com). Be sure to include the manual part number and the revision date listed on the front cover.

## Warranty

Thermo Scientific Laboratory Temperature Control Products have a warranty against defective parts and workmanship for 36 months from date of shipment. See back page of this manual for more details.

## Unpacking

Retain all cartons and packing material until the bath is operated and found to be in good condition. If the bath shows external or internal damage contact the transportation company and file a damage claim. Under ICC regulations, this is your responsibility.

If an immersion circulator was purchased without a bath, carefully remove the label and small piece of Styrofoam securing the float switch before using.

**Refrigerated baths should be left in an upright position for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲**

## Section 1 Safety

### Safety Warnings

Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your circulator. If you have any questions concerning the operation of your circulator or the information in this manual, please contact us. See inside cover for contact information.

**DANGER** indicates an imminently hazardous situation which, if not avoided, *will* result in death or serious injury.

**WARNING** indicates a potentially hazardous situation which, if not avoided, *could* result in death or serious injury.

**CAUTION** indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It also alerts against unsafe practices.

The lightning flash with arrow symbol, within an equilateral triangle, is intended to alert the user to the presence of non-insulated "dangerous voltage" within the circulator's enclosure. The voltage magnitude is significant enough to constitute a risk of electrical shock.

This label indicates the presence of hot surfaces.

This label indicates read the manual.

This label indicates a flammable refrigerant is used.

**Observe all warning labels. ▲**

**Never remove warning labels. ▲**

**Leave refrigerated baths in an upright position at room temperature (~25°C) for 24 hours before starting. This ensures the lubrication oil drains back into the compressor. ▲**

**The circuit protector located on the rear of the circulator is not intended to act as a disconnecting means. ▲**

**Operate the circulator using only the supplied line cord. The main power cord is used as the disconnecting device and must be easily accessible at all time. ▲**

Never operate the bath with the immersion circulator removed. ▲

Do not mount the immersion circulator backwards; the line cord could contact the reservoir fluid. Ensure the electrical cords do not contact any of the plumbing connections or tubing. ▲

Never place the circulator in a location or atmosphere where excessive heat, moisture, or corrosive materials are present. ▲

Ensure the tubing you select meets your maximum temperature and pressure requirements. ▲

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Never operate the circulator without fluid in the reservoir. ▲

Oil-based fluids expand when heated. Fill the reservoir only to the minimum level. ▲

Use only the approved fluids listed in this manual. Using other fluids voids the warranty. ▲

Other than water, before using any approved fluid, or when performing maintenance where contact with the fluid is likely, refer to the SDS and EC Safety Data Sheet for handling precautions. ▲

When using ethylene glycol and water, check the fluid concentration and pH on a regular basis. Changes in concentration and pH can impact system performance. ▲

Do not clean the circulator with solvents, use a soft cloth and water. ▲

Drain the reservoir before the bath is transported, and/or stored. ▲

Always turn the circulator off and disconnect the supply voltage from its power source before moving or before performing any service or maintenance procedures. ▲

Transport the circulator with care. Sudden jolts or drops can damage its components. ▲

Refer service and repairs to a qualified technician. ▲

The user is responsible for decontamination if hazardous materials are spilled. Consult the manufacturer regarding decontamination and or cleaning agents compatibility. ▲

Performance of installation, operation, or maintenance procedures other than those described in this manual may result in a hazardous situation and voids the manufacturer's warranty. ▲

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Many refrigerants which may be undetectable by human senses are heavier than air and will replace the oxygen in an enclosed area causing loss of consciousness. Contact with leaking refrigerant will cause skin burns. Refer to the bath's nameplate for the refrigerant used and then the refrigerant's SDS for additional information. ▲

Ensure, that no toxic gases can be generated by the fluid. Flammable gases can build up over the fluid during usage. ▲

Ensure the fluid is at a safe temperature (~40°C) before handling or draining. ▲

Never operate damaged or leaking equipment, or with any damaged cords. ▲

Never operate the circulator or add fluid to the reservoir with panels removed. ▲

## Personal Protective Equipment

Unless required by the manufacturer of the fluid used, there are no special personal protective equipment requirements needed to perform normal operation. Consult their Safety Data Sheet and/or EC Safety Data Sheet.. We do recommend wearing eye protection and gloves.

## Training

The user must review and understand all the sections in this manual before operating the bath.

## Element Assessment

### Chemical Used —

For refrigerated baths only:

R134A Refrigerant - CFC Free 3.5 ounces      Lubricating Oil:

R290 Refrigerant (propane), 4.1 ounces      POE – 243ml

Location: Refrigeration system      Location: Compressor

### Hazard Communication — None

**Electrical** — All electrical energy sources for the circulator/bath are provided by the user.

**Hazardous Energy Isolation** — The user is responsible to ensure that the circuit providing power has an energy isolation device.

**Electrical** Remove power from the circulator by turning the circuit protector off and disconnecting the line cord, see Section 3

### **Electrical Requirements.**

**Fluid** Before using any fluid refer to the manufacturer's SDS and EC Safety Data Sheet for handling precautions.

**Refrigeration** The bath's refrigeration system is sealed and requires no lockout/tagout.

Remove any isolation devices before starting the circulator.

**Chemicals** — Prior to decontamination procedures should be developed and employed to minimize or mitigate opportunities for exposure to chemicals present in or around the circulating bath.

### Radioactive Material — None

### Magnetic Fields — None

**Sampling and Analysis for External Contaminants from Other Sources or from Spill Events** — It is the responsibility of the user to identify the presence of all external chemicals before the decontamination. This should be accomplished by collecting samples and analyzing them for any suspected chemicals used or stored in areas adjacent to the bath.

**Oxygen Deficiency** — Many refrigerants which may be undetectable by human senses are heavier than air and will replace the oxygen in an enclosed area causing loss of consciousness. Refer to the bath's nameplate and the manufacturer's most current SDS for additional information. Test the atmosphere in a confined space area or area where gases are being used to purge the circulating bath. If necessary, use an air supplied respirator.

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
## **Fluorinated Greenhouse Gases**

The following information is included to comply with REGULATION (EU) No 517/2014 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 April 2014 on fluorinated greenhouse gases:

**This product contains fluorinated greenhouse gases in a hermetically sealed system.**

**If a leak in the sealed system is detected, the operator shall repair without undue delay.**

Refer to the F-Gas Declaration of Conformity for additional information.

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**Declaration of conformity with Article 14 of Regulation (EU) No 517/2014 of the European Parliament and of the Council**

We, Thermo Fisher Scientific (Asheville) LLC, F-Gas Portal Registration Number 23643, declare under our sole responsibility that when placing on the market pre-charged equipment, which we import to or manufacture in the Union, the hydrofluorocarbons contained in that equipment are accounted for within the quota system referred to in Chapter IV of Regulation (EU) No 517/2014 as:

We hold authorisation(s) issued in accordance with Article 18(2) of Regulation (EU) No 517/2014 and registered in the registry referred to in Article 17 of that Regulation, at the time of release for free circulation to use the quota of a producer or importer of hydrofluorocarbons subject to Article 15 of Regulation (EU) No 517/2014 that cover(s) the quantity of hydrofluorocarbons contained in the equipment.

The hydrofluorocarbons contained in the equipment have been placed on the market in the Union, subsequently exported and charged into the equipment outside the Union, and the undertaking that placed the hydrofluorocarbons on the market made a declaration stating that the quantity of hydrofluorocarbons has been or will be reported as placed on the market in the Union and that it has not been and will not be reported as direct supply for export in the meaning of Article 15(2)(c) of Regulation (EU) No 517/2014 pursuant to Article 19 of Regulation (EU) No 517/2014 and Section 5C of the Annex to Commission Implementing Regulation (EU) No 1191/2014.

January 1, 2017

\_\_\_\_\_  
Date

Mark Pearson, Director Global Regulatory Affairs

\_\_\_\_\_  
[name and position of legal representative]

\_\_\_\_\_  
[signature of legal representative]

Manufacturer:

Thermo Fisher Scientific (Asheville) LLC  
275 Aiken Road  
Asheville, NC 28804  
U.S.A.

EU Only Representative:

Thermo Electron LED GmbH  
Robert-Bosch-Strasse 1  
D-63505 Langenselbold  
Germany  
**VAT ID Number: DE 812 403 137**

## Section 2 General Information

### Description and Intended Use

The Thermo Scientific ADVANCED Series of Heated Immersion Circulators are used with refrigerated and heated baths. All immersion circulators can pump to an external system. All circulators have a digital display and easy-to-use touch pad, five programmable setpoint temperatures, acoustic and optical alarms, and offer adjustable high temperature protection.

The circulator is designed for use in a clean laboratory environment.

### ARCTIC Refrigerated/Heated Bath Circulator Specifications

<b>Stainless Steel Refrigerated/Heated Bath Circulators</b>			
	<b>A10</b>	<b>A10B</b>	<b>A25</b>
<b>AC150 Temperature Range</b> °C °F	-10 to 100 14 to 212	-10 to 100 14 to 212	-25 to 150 -13 to 302
<b>AC200 Temperature Range</b> °C °F	-10 to 100 14 to 212	-10 to 100 14 to 212	-25 to 200 -13 to 392
<b>Bath Volume</b> liters gallons	4 - 6 1.1 - 1.6	17 - 30 4.5 - 7.9	7 - 12 1.8 - 3.2
<b>Cooling Capacity</b> watts	240	250	500
<b>Refrigerant</b>	R134a	R134a	R134a
<b>Overall Dimensions</b> (H x W x D)* mm inches	670 x 220 x 414 26.4 x 8.7 x 16.3	509 x 429 x 913 20.0 x 16.9 x 35.9	749 x 273 x 483 29.5 x 10.7 x 19.0
<b>Work Area Dimensions</b> (D x W x L)* mm inches	150 x 137 x 124 5.9 x 5.4 x 4.9	200 x 297 x 365 7.9 x 11.7 x 13.4	200 x 173 x 184 8.0 x 6.8 x 7.2
<b>Net Weight</b> kg/lb	27.5/60.6	44.5/97.9	36.1/79.5
<b>Electrical Requirements**</b> (Voltage ±10%)	115 V/60 Hz		

\*See page 2-52. Add ~26 mm (1 inch) to D for drain fitting.

\*\*See Section 3 for additional information.

- Thermo Fisher Scientific reserves the right to change specifications without notice.



## ARCTIC Refrigerated/Heated Bath Circulator Specifications

Stainless Steel Refrigerated/Heated Bath Circulators		
	A25B	A45HC
<b>AC150 Temperature Range</b> °C °F	-25 to 150 -13 to 302	
<b>AC200 Temperature Range</b> °C °F	-25 to 200 -13 to 392	-45 to 200 -49 to 392
<b>Bath Volume</b> liters gallons	13 - 21 3.4 - 5.5	7 - 12 1.8 - 3.2
<b>Cooling Capacity</b> watts	500	900
<b>Refrigerant</b>	R134a	R290
<b>Overall Dimensions</b> (H x W x D)* mm inches	778 x 324 x 541 30.6 x 12.8 x 21.3	787 x 385 x 580 31.0 x 15.2 x 22.8
<b>Work Area Dimensions</b> (D x W x L)* mm inches	233 x 224 x 244 9.2 x 8.8 x 9.6	200 x 173 x 184 8.0 x 6.8 x 7.2
<b>Net Weight</b> kg/lb	42.3/93.1	52/115
<b>Electrical Requirements**</b> (Voltage ±10%)	115 V/60 Hz	




\*Add ~26 mm (1 inch) to D for drain fitting.

\*\*See Section 3 for additional information.

- Thermo Fisher Scientific reserves the right to change specifications without notice.

### Cooling Capacity





Specifications obtained at sea level using water (above +5°C to +90°C) or a fluid with a specific heat of 2.3 kJ/kg-K or 0.55 Btu/lb-F (less than 5°C) as the recirculating fluid at a +20°C ambient condition, at nominal operating voltage. Other fluids, process temperatures, ambient temperatures, altitude or operating voltage will affect performance. The bath was operating with the lid closed and no heat load lid or any external circulation.

Minimum temperatures are only achieved with the above conditions. If your conditions vary from above, minimum temperature will most likely not be achieved.

Specifications are for reference only and are subject to change.

## Time to Temperature

Specifications obtained at sea level using water (above +5°C to +90°C) or a fluid with a specific heat of 2.3 kJ/kg-K or 0.55 Btu/lb-F (less than 5°C) as the recirculating fluid at a +20°C ambient condition, at nominal operating voltage. Other fluids, process temperatures, ambient temperatures, altitude or operating voltage will affect performance. Pump specifications are nominal values of ±10%. Specifications are for reference only and are subject to change

## SAHARA Heated Bath Circulator Specifications

Stainless Steel Bath Circulators		
	S7	S13
<b>AC150 Temperature Range</b> °C* °F*	Ambient +13 to 150 Ambient +23 to 302	Ambient +13 to 150 Ambient +23 to 302
<b>AC200 Temperature Range</b> °C* °F*	Ambient +13 to 200 Ambient +23 to 392	Ambient +13 to 200 Ambient +23 to 392
<b>Bath Volume</b> liters gallons	4 - 8 1.1 - 2.1	7 - 12 1.8 - 3.2
<b>Overall Dimensions</b> (H x W x D)** mm inches	494 x 235 x 428 19.5 x 9.2 x 16.7	494 x 321 x 428 19.5 x 12.6 x 16.7
<b>Work Area Dimensions</b> (D x W x L)** mm inches	200 x 154 x 112 7.9 x 6.1 x 4.4	200 x 293 x 112 7.9 x 9.4 x 4.4
<b>Net Weight</b> kg/lb	10.6/23.4	12.3/27.0
<b>Electrical Requirements***</b> (Voltage ±10%)	115 V/60 Hz	

Stainless Steel Bath Circulators			
	S21	S45	S49
<b>AC150 Temperature Range</b> °C* °F*	Ambient +13 to 150 Ambient +23 to 302	Ambient +13 to 150 Ambient +23 to 302	Ambient +13 to 150 Ambient +23 to 302
<b>AC200 Temperature Range</b> °C* °F*	Ambient +13 to 200 Ambient +23 to 392	Ambient +13 to 150 Ambient +23 to 302	Ambient +13 to 200 Ambient +23 to 392
<b>Bath Volume</b> liters gallons	7 - 19 1.8 - 5.0	30 - 41 7.9 - 10.8	29 - 53 7.7 - 14.0
<b>Overall Dimensions**</b> (H x W x D) mm inches	447 x 381 x 628 17.6 x 15.0 x 24.7	594 x 381 x 628 23.4 x 15.0 x 24.7	494 x 579 x 746 19.5 x 22.8 x 29.4
<b>Work Area Dimensions</b> (D x W x L) mm inches	150 x 297 x 312 5.9 x 11.7 x 12.3	300 x 298 x 312 11.8 x 11.7 x 12.3	200 x 498 x 430 7.9 x 19.6 x 16.9
<b>Net Weight</b> kg/lb	14.2/31.2	20.3/44.7	24.3/53.4
<b>Electrical Requirements***</b> (Voltage ±10%)	115 V/60 Hz		



\*The lowest usable temperature may be 13°C (23°F) warmer than room temperature, lower temperature ranges require supplemental cooling.

\*\*See page 2-2. Add ~26 mm (1 inch) to D for drain fitting.

\*\*\*See Section 3 for additional information.

- Thermo Fisher Scientific reserves the right to change specifications without notice.

**Time to Temperature**

Specifications obtained at sea level using water (above +5°C to +90°C) or a fluid with a specific heat of 2.3 kJ/kg-K or 0.55 Btu/lb-F (less than 5°C) as the recirculating fluid at a +20°C ambient condition, at nominal operating voltage. Other fluids, process temperatures, ambient temperatures, altitude or operating voltage will affect performance. Pump specifications are nominal values of ±10%. Specifications are for reference only and are subject to change. Heat-up rates for the 100V baths will take approximately 25% longer than the 115V.

**Wetted Materials****ADVANCED Immersion  
Circulator**

Viton  
EPDM  
Ceramic  
Ultem  
Vectra  
Stainless Steel

**Stainless Steel Baths/Circulators**

Stainless Steel 316  
Stainless Steel 304  
EPDM  
Ryton (drain fitting)  
Zotek-N (cover seal)

## Section 3 Installation

### Ambient Conditions

Ambient Temperature Range	5°C to 40°C (41°F to 104°F)
Maximum Relative Humidity	80% for temperatures up to 31°C (88°F)*
Operating Altitude	Sea Level to 2000 meters (6560 feet)
Overvoltage Category	II
Pollution Degree	2
Degree of Protection	IP 20

\*decreasing linearly to 50% relative humidity at 40°C (104°F)


Never place the circulator in a location where excessive heat, moisture, inadequate ventilation, or corrosive materials are present. ▲

Carefully install the circulator to ensure it does not fall into the bath or that its line cord does not make contact with the bath contents. ▲

### Ventilation

All circulators except A45HC can operate with 0 clearance on two exhaust sides as long as the third exhaust side has unrestricted air flow. A45HC must have at least 25mm [1"] of clearance on both sides to comply with UL leak mitigation safety standards. Blocked ventilation increases the circulator's temperature, reduce its cooling capacity and, on refrigerated baths, eventually lead to premature compressor failure.

Ventilation Options


## Bath Circulator

Circulators are designed for continuous operation; and for indoor use.

The equipment normally ships with the immersion circulator mounted facing the reservoir. You may change the position  $\pm 90^\circ$  by removing the thumb screws, no tools are required.

**Do not mount the circulator backwards; the line cord could contact the reservoir fluid. ▲**

**Never place the circulator in a location where excessive heat, moisture, inadequate ventilation, or corrosive materials are present. ▲**

**Leave refrigerated baths in an upright position at room temperature ( $\sim 25^\circ\text{C}$ ) for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲**

**A45HC units utilize a flammable refrigerant, risk of fire or explosion, do not use sharp objects to clean evaporator coil in the tank. ▲**

## Electrical Requirements

The circulator's construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided. ▲

For Sahara baths, and stand-alone immersion circulators, use the supplied nut and washer to secure the grounding strap to the top of the metal tank. See figure below. ▲

The circulator is intended for use on a dedicated outlet. All circulators are equipped with automatic thermally-triggered 20 Amp circuit protector.

The circuit protector is designed to protect the circulator.

**Note** If the circuit protector activates allow the circulators to cool before resetting. Restart the circulator. Contact us if it activates again. ▲

The circulator's power cord is used as the disconnecting device, it must be easily accessible at all times. ▲


Refer to the bath nameplate on the rear, upper-left-hand corner of the bath for specific electrical requirements. Voltage deviations of  $\pm 10\%$  are permissible. The outlet must be rated as suitable for the total power consumption of the circulator, see next page.

## Grounding Strap and Nut

Sahara baths, and stand-alone circulators, require immersion circulator grounding.


For Sahara baths, and stand-alone immersion circulators, use the supplied nut and washer to secure the grounding strap to the top of the metal tank. ▲


The following power options are available:

<b>...Bath</b>	<b>Volts<sup>1</sup>/Hertz/Phase</b>	<b>Amps<sup>2</sup></b>	<b>Total Wattage</b>	<b>Plug Type</b>
<b>A10</b>	115/60/1	11.5	1165	N5-15
<b>A25/A25B</b>	115/60/1	11.7	1185	N5-15
<b>A10B</b>	115/60/1	11.5	1165	N5-15
<b>A45HC</b>	115/60/1	14.4	1660	N5-20
<b>All Heated Baths/Circulators</b>	115/60/1	11.3	1300	N5-15

1. Volts  $\pm$  10%
2. Maximum amp draw

20 Amp Outlet  
(16 Amp)


15 Amp Outlet  
(12 Amp)

**For refrigerated baths:**

Ensure all communication and electrical connections are made prior to starting the circulator. ▲

- Install the supplied RJ45 shielded cable between the immersion circulator and the bath RJ45 connectors (similar to Ethernet). **This is required for proper operation.**


- Install the power cord from the connector on the rear of the controller, A, to the connector on the rear of the refrigerated bath, B.
- Connect the bath's power cord, C, to a grounded power outlet.


For refrigerated baths, never connect controller power inlet, A, to a power outlet. Never connect power outlet, B, to anything but an immersion circulator. ▲

Ensure the electrical cords do not come in contact with any of the plumbing connections or tubing. ▲


**For non-refrigerated baths:**

- Connect the bath's power cord, C, to a grounded power outlet.


**Ensure the electrical cords do not come in contact with any of the plumbing connections or tubing. ▲**

## External Circulation

The plumbing connections for external circulation are located on the rear of the immersion circulator.  is the return flow from the external application.  is the outlet flow to the external application (supply side). The connections are 16 mm O.D. Remove the union nuts and plates to install the 8 mm or 12 mm hose barbs and clamps supplied with the circulator.

To prevent damage to the circulator's plumbing, use a 19 mm backing wrench when removing/installing the connections. ▲

## Remote Temperature Sensor

The remote temperature sensor on the rear of the immersion circulator requires a 4-pin connector that must mate to a LEMO # ECP.1S.304.CLL. The immersion circulator uses a 3 wire sensor, but a 4 wire sensor can be used (pins 3 and 4 are interconnected in the control head). The pin-out is:

Pin 1 and 2 = Pt100 +      Pin 3 and 4 = Pt100 -

See Section 4 for instructions to enable the remote sensor.

**Note** Bath circulator accessories are not part of the UL Listings. ▲

## Multifunction Port

AC 200s have a multifunction port on the rear of the immersion circulator that is used in conjunction with an *optional* accessory box to select and activate additional features. See Section 4 for additional information.

**Note** Bath circulator accessories are not part of the UL Listings. ▲

## USB Port

AC 200 only. If your computer does not automatically recognize the USB driver, troubleshooting instructions are provided in Section 6.

## Tubing Requirements

**Ensure none of the tubing comes in contact with the power cord.** ▲

Tubing is normally used to connect the pump to an external application.

**Note** The maximum allowable length of tube depends largely on the size, form and material of the external vessel. The length of tube and its diameter, combined with the circulating capacity, have a large affect on the temperature stability. Whenever possible, use a wider tube diameter and place the application as close as possible to the circulator. ▲

**Extreme operating temperatures will lead to extreme temperatures on the tube surface, this is even more critical with metal nozzles.** ▲

- the required tube material depends on the heat transfer liquid used
- tubes must not be folded or bent
- after prolonged use, tubes may become brittle or they may get very soft, check them on a regular basis and replace if necessary
- secure all tube connections using clamps

When using the internal bath only, the plumbing connections can be closed with the supplied plate and union nuts.

## Tubing

Tubing for Thermo Scientific temperature control systems is optional. Please select the proper tubing from the table shown in Section 5.

**Ensure the tubing you select meets your maximum temperature and pressure requirements. ▲**

### Plastic and rubber tubing

If other plastic and rubber tubes are used, ensure that the tubes selected are fully suitable for the particular application, i.e., that they will not split, crack or become disengaged from their connections.

Connect the tubing using the supplied tube fittings for 8 or 12 mm i.d. They are attached to the plumbing connections with a supplied coupling nut.

We highly recommend using foam rubber insulation on the tubing and the fittings.

### Metal tubing

Thermo Scientific metal tubing (stainless steel insulated) offers a particularly high degree of safety and is suitable for both low and high temperatures/liquids.

The metal tubing is attached directly to the plumbing connections, gaskets are not required.

**Do not subject tubing to mechanical strain and ensure any specified bend radius is not exceeded. ▲**

Tubing is available in lengths of 0.5, 1.0 and 1.5 meters. Couplings for connecting tubes are also available.

The smallest opening inside the metal tubes is 10 mm. The metal tubing is provided with coupling nuts (M16 x 1, DIN 12 879, part 2) at either end.

## Approved Fluids

The user is always responsible for the fluid used. Never use corrosive fluids with this circulator. ▲

Never use 100% glycol. ▲

Handle and dispose of liquids, other than water, in accordance with the fluid manufacturers specification and/or the SDS. ▲

Adjust the circulator's software to the fluid used, see Section 4. ▲

When using water above 80°C closely monitor the fluid level, frequent top-offs will be required. It also creates steam. ▲

Water/glycol mixtures require top-offs with pure water, otherwise the percentage of glycol will increase resulting in high viscosity and poor performance. ▲

The fluid information is only a guide since specifications may change.

Thermo Fisher Scientific takes no responsibility for damages caused by the selection of an unapproved bath fluid.

Unapproved bath fluids are fluids which:

- are very highly viscous (much higher than 30 mPas at the respective working temperature)
- have corrosive characteristics or
- tend to break down at high temperatures

For fluid selection consider application requirements, operating temperature range, material compatibility, safety concerns, and environmental issues.

### Chlorine

Short term usage of tap water may not cause any adverse affects on the circulator or your application, but in the long term problems may arise. To help alleviate these problems Thermo Fisher Scientific recommends the use of chlorine.

The duration of time that chlorine remains in solution depends on factors such as water temperature, pH and availability of direct sunlight. We recommend maintaining chlorine levels at proper levels using chlorine test strips, generally 1 to 5 ppm is adequate.

For best results, maintain the pH of the fluid between 6.5 and 7.5. Do not add additional chlorine without first determining the concentration ratio

that already exists in the fluid supply. Corrosion and degradation of the circulation components can result from concentration ratios that are too high. Contact our customer support for additional information.

**5°C to 95°C — Distilled Water or Deionized Water** (up to 3 MΩ-cm)  
Normal tap water leads to calcareous deposits necessitating frequent circulator decalcification, see table on next page.

Calcium tends to deposit itself on the heating element. The heating capacity is reduced and service life shortened.

**-30°C to 80°C — Water with Glycol**

Below 5°C water has to be mixed with a glycol. The amount of glycol added should cover a temperature range 5°C lower than the operating temperature of the particular application. This prevents the water/glycol from gelling (freezing) near the evaporating coil.

Excess glycol deteriorates the temperature accuracy due to its high viscosity.

**-40°C to 200°C — SIL180:**

SIL180 is suitable for covering nearly the entire range with just one liquid, especially when used with the cooling circulators. See next page for additional information.

Unfortunately SIL180 has a wetting tendency necessitating the occasional cleaning of the bath cover.

**other temperatures:**

Thermo Fisher Scientific offers a range of heat transfer fluids for these temperature control applications.

**SYNTH 60 and SYNTH 260:**

Synthetic thermal liquid with a medium life span (several months) and little smell annoyance. See next page for additional information.

**SIL 100, SIL 180, SIL 300:**

Silicone oil with a very long life span (over 1 year) and negligible smell. See next page for additional information.

Thermo Fisher heat transfer fluids are supplied with an EC Safety Data Sheet.

## Additional Fluid Precautions

Ensure, when selecting the heat transfer fluid, that no toxic gases are generated. Flammable gases can build up over the fluid during usage. ▲

Ensure the circulator's High Temperature Cutout (HTC) is set to 25°C below the fire point of the selected fluid. ▲

The highest working temperature must be limited to 25°C below the fire point of the bath fluid, as defined by UL/CSA/EN/IEC 61010-1 and 61010-2-010. ▲

When working with fluids other than water:

- Do not use any fluid until you have read and understood the label and the Safety Data Sheet (SDS).
- Do not blend any fluids.
- Ensure any fluid residue or any other material is thoroughly removed before filling the bath with a different fluid.
- Always wear protective clothing, especially a face shield and gloves.
- Avoid spattering on the circulator's components, always *slowly* add fluid. When adding fluid, point the opening of a container away from yourself.
- Use fume hoods.
- Do not allow any ignition sources in the vicinity.

### Fire point

Flammable thermal liquids can ignite when a specified temperature is surpassed. The bath liquid is limited to a temperature level 25°C below the fire point as defined by UL/CSA/EN/IEC 61010-1 and 61010-2-010.

### Viscosity

For optimum temperature accuracy, it is important that heat transfer liquids have a low viscosity.

#### ■ Working temperature range

This is the recommended long-term operating range. The maximum viscosity is approximately 5 mPas.

#### □ Operating temperature range


Long-term operation is recommended only under certain conditions. The viscosity may rise to a maximum of 30 mPas. The pump capacity will not match specifications.

#### ■ Heating-up range

Long-term operation is not recommended, the pump motor's excess temperature protection may switch off the pump.

### Range of Application

	Sil 100	Sil 180	Sil 300	Synth 60	Synth 260
Fire Point °C	>100	>225	>325	70	275
Flash Point °C	57	170	300	59	260
Viscosity at 20°C (mPas)	3	11	200	2	140
Density at 20°C (kg/dm³)	0.89	0.93	1.08	0.76	1.03
Specific heat capacity (kJ/Kg*K)	1.67	1.51	1.56	2.10	2.00

Order Number 10 liter container	999-0202	999-0204	999-0206	999-0210	999-0214
Order Number 5 liter container	999-0201	999-0203	999-0205	999-0209	999-0213
Color	trans-parent, colorless	trans-parent, colorless	trans-parent, colorless	trans-parent, colorless	trans-parent, yellow
Reacts with	Silicone	Silicone	Silicone	Rubber Silicone	Copper Light metals Bronze

EC-Safety Data Sheets will be delivered together with each container of liquid.

## Water Quality and Standards

Process Fluid	Permissible (PPM)	Desirable (PPM)
<b>Microbiologicals</b>		
(algae, bacteria, fungi)	0	0
<b>Inorganic Chemicals</b>		
Calcium	<25	<0.6
Chloride	<25	<10
Copper	<1.3	<1.0
	0.020 ppm if fluid in contact with aluminum	
Iron	<0.3	<0.1
Lead	<0.015	0
Magnesium	<12	<0.1
Manganese	<0.05	<0.03
Nitrates\Nitrites	<10 as N	0
Potassium	<20	<0.3
Silicate	<25	<1.0
Sodium	<20	<0.3
Sulfate	<25	<1
Hardness	<17	<0.05
Total Dissolved Solids	<50	<10
<b>Other Parameters</b>		
pH	6.5-8.5	7-8
Resistivity	0.01*	0.05-0.1*

\* MΩ-cm (compensated to 25°C)

Unfavorably high total ionized solids (TIS) can accelerate the rate of galvanic corrosion. These contaminants can function as electrolytes which increase the potential for galvanic cell corrosion and lead to localized corrosion such as pitting. Eventually, the pitting will become so extensive that refrigerant will leak into the water reservoir.

As an example, raw water in the United States averages 171 ppm (of NaCl). The recommended level for use in a water system is between 0.5 to 5.0 ppm (of NaCl).

Recommendation: Initially fill the tank with distilled or deionized water. Do not use untreated tap water as the total ionized solids level may be too high. This will reduce the electrolytic potential of the water and prevent or reduce the galvanic corrosion observed.

## Filling Requirements

Ensure the reservoir drain port on the front of the bath is *closed* and that all plumbing connections are secure. Also ensure any residue is thoroughly removed before filling the bath.

**Before using any fluid refer to the manufacturer's SDS and EC safety data sheets for handling precautions. ▲**

To avoid spilling, place your containers into the bath before filling.

With a low level WARNING the circulator continues to run, with a FAULT the circulator shuts down the refrigeration, pump and heater, see Section 6. The low level warning is at approximately 43 mm (1 3/4") below the top, the low level fault is at approximately 57 mm (2 1/4").

**Oil-based fluids expand when heated. Fill the reservoir only to the minimum level. Monitor the fluid level whenever heating the fluid. ▲**

When pumping to an external system, keep extra fluid on hand to maintain the proper level in the circulating lines and the external system.

## Draining

**Before draining any fluid refer to the manufacturer's SDS and EC safety data sheets for handling precautions. ▲**


**Ensure the fluid is at a safe handling temperature, ~40°C. Wear protective clothing and gloves. ▲**

- place a suitable vessel underneath the drain. If desired, attach an 8 mm id tube on the drain.
- *slowly* turn the drain plug until flow is observed.

**Turning the drain cap more than 1 1/2 turns will result in the drain cap and fitting coming off the bath. ▲**

In this case screw the drain fitting back into the bath. Attaching the cap onto the fitting aids in installation. If required, contact us for additional information.


Installed Drain Fitting with Cap Removed

## Section 4 Operation

### ADVANCED Heated Immersion Circulator

The Thermo Scientific ADVANCED Series of Heated Immersion Circulators have a digital display and easy-to-use touch pad, five programmable setpoint temperatures, acoustic and optical alarms. All circulators offer adjustable high temperature protection.

This label indicates read the instruction manual before starting the circulator.

Use this button to place the circulator in and out of standby, see page 4-3 for more details. The blue LED illuminates when standby is enabled.

Use these navigation arrows to move through the circulator displays and to adjust values.


Pressing this button once to make changes on the immersion circulator's display screen. In most cases, pressing it again is required to save the change.

Use this button to cancel any changes and to return the immersion circulator to its previous display. Canceling a change can only be made before the change is saved. In some cases, it is also used to save changes.

**Note** Holding this button for five seconds resets the display contrast to the default level and also brings up the language menu to change, if needed, the displayed language. See **Settings-Display Options** in this Section. ▲

Used for adjusting and resetting the High Temperature Cutout. Details are explained in this Section.

## Setup

Leave refrigerated baths in an upright position at room temperature ( $\sim 25^{\circ}\text{C}$ ) for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲

Before starting, double check all USB (optional), electrical and plumbing connections. ▲

Ensure the circulator's HTC is set to  $25^{\circ}\text{C}$  below the fire point of the selected fluid. ▲

## Start Up


Do not run the circulator until fluid is added to the reservoir. Have extra fluid on hand. If the circulator does not start refer to Section 6 Troubleshooting.

- On refrigerated baths, place the circuit protector located on the rear to the **I** position.
- Place the circuit protector located on the back of the immersion circulator to the **I** position.

After a slight delay the blue LED on the front panel illuminates.

- Press , the Start Display appears and the blue LED goes out.
- Ensure the start symbol has a highlight box around it, if not use the arrows to navigate to the symbol.

- Press . The circulator starts and the start symbol turns into a stop symbol (■).

**Note** The pump starts immediately but refrigerated baths take up to 30 seconds before the compressor starts. ▲


**Note** After start up, check all the plumbing connections for leaks. ▲

The **SP1** and **Menu** portions on the top of the display are used to view and/or change the circulator's settings. They are explained in detail later in this Section.

**In** indicates the circulator is using its internal sensor for temperature control. **Ex** indicates the circulator is using an external sensor for temperature control.

## Status Display

If desired, press  to toggle between the Start/Status Displays.

**Note** If no operator inputs are being made, the circulator automatically switches to the Status Display after 60 seconds. If desired, change the time or disable this feature using the **Display Options** Menu. ▲


## Stand By Mode

Press , the display goes blank and the circulator enters the stand by mode. The blue LED illuminates.

## Stopping the Circulator

Ensure the stop symbol is highlighted, if not use the arrows to navigate to the symbol.

Press . The circulator stops and the stop symbol turns into a start symbol ().

## Power Down

Press , the display goes blank and the circulator enters the stand by mode. The blue LED illuminates.

## Shut Down

Place the circuit protector on the back of the circulator to the  position. The blue LED extinguishes.

On refrigerated baths, place the circuit protector on the rear of the bath to the  position.

Using any other means to shut down a refrigerated bath can reduce the life of the compressor. ▲

Always turn the circulator off and disconnect it from its supply voltage before moving it. ▲

The circuit protector located on the rear of the components is not intended to act as a disconnecting means. ▲

## Restarting

**Note** When quickly restarting, the compressor may take up to 10 minutes before it starts to operate. ▲


## Changing the Setpoint

**Note** You can not adjust the setpoint closer than 0.1°C to either of the fluid's system limits, see Fluids Type in this Section, or beyond the circulator's temperature range. ▲

The setpoint can be changed with the circulator running or not.


The Setpoint is the desired process fluid temperature. The circulator can store up to five setpoints, **SP1** through **SP5**. The procedure for changing stored setpoint values is discussed later in this Section.

Use the navigation arrows and move to the setpoint window and then press  to highlight it and the setpoint value as shown below.


Use the up and down navigation arrows to bring up the desired setpoint, **SP1** through **SP5**, and then press .

The Setpoint Value Window now indicates the corresponding setpoint's stored value.

If desired, you can change the displayed setpoint value by using the navigation arrows to highlight the Setpoint Value Window and then pressing . The right-most digit will have a cursor beneath it.

Use the left and right arrows to move the cursor to the desired digit and then use the up and down arrows to change the value. Once all the desired changes are made, press  to save the change.

**Note** Using this procedure also changes the setpoint's stored value. ▲


## Menu Displays

The circulator uses menus to view/change the settings.


**Note** The circulator does not need to be running to view/change these settings. ▲

For all **Menu** displays, once  is pressed to change a display, you can press  to return to the previous screen.


1. Use the arrow buttons to highlight **Menu** and the circulator brings up the Main Menu Display.
2. Use the up and down arrow to highlight the desired setting and then press  to bring up additional submenus.

Main Menu Display

See page 4-8.



See page 4-17.

Since the circulator can only display five lines of text at a time, keep pressing the down arrow to view any additional options.

## Menu

The **Menu** line, at the bottom of all the submenu displays, is another way to return the circulator back to the Start Display.


1. From any submenu display, use the down arrow button to highlight **Menu**.
2. Press  to return to the Start Display.

Start Display

Pressing  from the Menu line returns you to the previous screen.



# Menu Tree

Since the circulator can only display five lines of text at a time, keep pressing the down arrow to view any additional options.

**Settings - Application Settings** is used to view/adjust the circulator's five Setpoints (**SP**) and Real Temperature Adjustments (**RTA**), enable/disable the alarms, change the fluid type, set the pump speed, configure the interface (optional), set the clock, turn the timer on or off, and turn auto restart and energy savings on or off.



1. With **Application Settings** highlighted press  to view:
2. Scroll down for additional options.

3. With **Setpoints** highlighted, press  to display the list. Use the up/down arrows to highlight the desired **SP**. (Use the down arrow to display **SP5**.)

4. Press .

The Setpoint and RTA are changed using the same procedure. With the desired setpoint highlighted press  to display the submenu.

If this temperature on the Start/Status Displays does not accurately reflect the actual temperature in the reservoir, an RTA can be applied. The RTA can be set  $\pm 10^{\circ}\text{C}$  ( $\pm 18^{\circ}\text{F}$ ).

As an example, if the circulator's temperature is stabilized and displaying  $20^{\circ}\text{C}$  but a calibrated reference thermometer reads  $20.5^{\circ}\text{C}$ , set the RTA to  $-0.5^{\circ}\text{C}$ . After you enter a RTA value allow the temperature to stabilize before verifying the temperature in the bath. **Note** If display accuracy is required, we recommend repeating this procedure at various setpoint temperatures and on a regular basis. ▲

**Note** You cannot adjust the setpoint closer than  $0.1^{\circ}\text{C}$  to either of the fluid's system limits, see Fluids Type in this Section. ▲


5. With the desired line highlighted press .

The right-most digit will have a cursor beneath it. Use the left and right arrows to move the cursor to the desired digit and then use the up and down arrows to change the value. Once all the desired changes are made, press  to save the change or  to cancel it.


**35.0<sup>°</sup>C**

**Alarms** is used to view/adjust the high and low temperature alarm limits, to enable/disable the audible alarms and to configure the low level warning reaction. **Note** The alarm range limits depend on the fluid, circulator and bath used, see next page. ▲

1. With **Alarms** highlighted, press  to display:

2. With **Temperature Alarms** highlighted, press  to display:

3. Highlight the desired limit and press . Follow the same procedure used to change a setpoint.

If the **Fault** temperature is exceeded the circulator shuts down and, if enabled, the audible alarm sounds. If the **Warn** temperature is exceeded the circulator continues to run and, if enabled, the audible alarm sounds. In either case a message is displayed.

**High Fault** cannot be set below **High Warn**.


**High Warn** cannot be set below **Low Warn**.

**Low Fault** cannot be set above **High Warn**.

**Note** When changing the temperature alarms the current setpoint also changes if it falls outside the new limits. ▲

1. With **Audible Alarms** highlighted, press  to display the alarms. The factory presets are off.

Highlight the desired alarm and press  to toggle between enable and disable mode.


If **Faults** is enabled the alarm sounds when a fault occurs. If **Warnings** is enabled the alarm sounds when a warning occurs.

If **Prog. End** is enabled the circulator beeps twice at the end of each cycle and three times at the end of the program. If **Prog. Step** is enabled the circulator beeps once at the beginning of the program and once at the end of each step. Press  to return to the previous display.

1. With **Low Level** highlighted, press  to toggle the low level warning alarm on/off.

The factory preset is off.

Press  to return to the previous display.


**Note** The low level cutout is always enabled and if activated will generate a fault message. ▲

**Fluids Type** is used to identify the type of fluid used. The circulator uses the fluid type to automatically set certain operating parameters.

1. With **Fluid Type** highlighted, press  to display the list of approved fluids.

Highlight the desired fluid and then press  to select it.


2. With the desired fluid selected press  to return to the previous display.

**Note** The circulator's operating range is determined by the currently selected fluid. When a new fluid is selected the circulator, if necessary, automatically adjusts the temperature alarms and/or setpoint. ▲

**Pump Adjustment** is used to review/set the desired pump speed.

1. With **Pump Adjustment** highlighted, press  to display the speeds.

Highlight the desired speed and press  to select it.

	Fluid temperature alarm limits	
	High °C	Low °C
<b>AC150 circulators:</b>		
Water	+95	+5
Ethylene Glycol-Water	+100	-30
Propylene Glycol-Water	+100	-20
Other	+150	-40
SIL 100	+75	-40
SIL 180	+150	-40
SIL 200	+150	+20
SIL 300	+150	+80
SYNTH 60	+45	-40
SYNTH 260	+150	+45


<b>AC200 circulators:</b>		
Water	+95	+5
Ethylene Glycol-Water	+100	-30
Propylene Glycol-Water	+100	-20
Other	+200	-40
SIL 100	+75	-40
SIL 180	+200	-40
SIL 200	+200	20
SIL 300	+200	+80
SYNTH 60	+45	-40
SYNTH 260	+200	+45

**Note** The range is also limited by the bath temperature range, see Section 2. ▲

**Note** When using oil as a reservoir fluid, we recommend running the circulator at 95°C for 15 minutes to remove any moisture in the fluid. ▲


**Set Clock** is used to set the circulator's **Set Time/Date** (hr : min : sec) and date (year - month - day).

**Format Date** is only applied to the date sent out the serial interface **DD/MM/YYYY** or **MM/DD/YYYY**. The date displayed is fixed at **year - month - day**.

**On/Off Timer** is used to enable and set the circulator's timer.

1. With **On/Off Timer** highlighted, press  to display the on (I) and off (O) time as well as the enable box.

After setting the on and off times highlight **Enable** and press  to activate the timer.

**Remote Sensor** is used to enable the optional remote temperature sensor feature, see Section 3.

1. With optional **Remote Sensor** highlighted, press  to toggle between enable and disable.

**Auto Restart** is used to enable the auto restart feature. When enabled, the immersion circulator automatically restarts after a power failure or power interruption condition is restored. If a ramp was running when power failed, the ramp program resumes where it left off. **Note** Consider any possible risks before enabling this mode of operation. ▲

1. With **Auto Restart** highlighted, press  to toggle between enable and disable.

**Energy Saving** is used to enable the energy savings mode. The Energy Saving mode is primarily designed for applications running under a stable load. Enabling the mode saves energy by reducing the circulator's heater power and cooling requirements. This can result in substantial energy savings over the life of the circulator. The default setting is disable (enable for A45HC refrigerated circulators).

1. With **Energy Saving** highlighted, press  to toggle between enable and disable.

## Edit - Ramp Program (AC 200 only)

**Edit Ramp Program** is used to view/adjust the immersion circulator's program function.

Define your program as a series of setpoints with a known period of time interval between each. Each interval is one step of the program. Pay careful attention to the first part of your program. What conditions must exist at the beginning of your process? For example, at the starting setpoint you may need to program an initial period of constant temperature to allow for thermal stabilization.

**Note** Consider the circulator's limitations when designing programs. Temperature or time parameters which exceed the performance capabilities of the circulator will result in unsatisfactory operation. If reaching the ramp setpoint temperatures is important, you will have to operate the bath between the desired setpoints and note the duration before programming the ramp. ▲

It is possible to create a program calling for very rapid changes in temperature. Although the circulator may not be capable of producing such changes, it may be practical to program such steps as a way to cause the fastest possible temperature change.

1. With **Ramp Program** highlighted press  to display:

<b>Edit Ramp</b>	▲
<b>Edit Step</b>	
	▼
<b>Menu</b>	

2. With **Edit Ramp** highlighted press  to display:

<b>No of Steps</b>	<b>XX</b>	▲
<b>Variance</b>	<b>xxx.xx</b>	
<b>Cycles</b>	<b>xxx</b>	
<b>End State</b>	<b>Shut Down</b>	▼
<b>Menu</b>		

The ramp can have up to 30 **Steps**.

The **Variance** is used to set a temperature range, the program starts when the fluid temperature is within this range. For example, if the desired **Start Temp** is 25°C and the **Variance** is set to +5°C, the program automatically starts when the bath temperature is between 20°C to 30°C.

**Cycles** sets the number of times the entire ramp program is repeated after the last step is completed. For example, selecting 3 **Cycles** runs the entire ramp program a total of 3 times.

**End State** configures the circulator to either **Shut Down** or continue running (**Maintain**) when the program is over.

3. Once the **Edit Ramp** portion is complete press  and then highlight **Edit Step**.

<b>Edit Ramp</b>	^
<b>Edit Step</b>	
	v
<b>Menu</b>	

4. Use **Edit Step** to enter the parameters for each step. The maximum **Duration** for a step is 1000 minutes.

<b>Step #</b>	<b>XX</b>	^
<b>Start Temp</b>	<b>xxx.xx</b>	
<b>Stop Temp</b>	<b>xxx.xx</b>	
<b>Duration (min)</b>	<b>xxxxx</b>	v
<b>Menu</b>		

Other than step one, the circulator automatically sets the **Start Temp** with the **Stop Temp** from the previous step. This value cannot be changed.

The ramp program has an **Assured Soak** feature that can be enabled for each step independently. When enabled this feature pauses the ramp timer until the temperature reaches setpoint,  $\pm$  variance. This assures the temperature reaches setpoint,  $\pm$  variance, before the ramp program continues to the next step.

5. Scroll down to view the **Assured Soak** feature. Press  to enable/disable.

<b>Start Temp</b>	<b>xxx.xx</b>	^
<b>Stop Temp</b>	<b>xxx.xx</b>	
<b>Duration (min)</b>	<b>xxxxx</b>	
<input type="checkbox"/> <b>Assured Soak</b>		v
<b>Menu</b>		

6. After all the desired steps are built, keep pressing  until the Start Display appears.


<b>SP1</b>	◇	<b>Menu</b>
<b>20.00°C</b>		
		(In)
<b>24.29 °C</b>		

You can enable an alarm to sound when each step and/or the program is complete, see **Settings - Basic Settings** in this section.



## Running a Ramp Program


1. Use the arrows to highlight **SP** window and press .

2. Press the up or down arrow until the window displays **Ramp**.

3. Press  and the circulator displays the start  and pause/resume  program symbols. Highlight the start symbol .

4. With the start symbol highlighted, press  to start the program. The start symbol changes into a highlighted stop program symbol .


**Note** With the stop symbol  highlighted, pressing  shuts down the circulator. ▲

The program does not start until the process fluid temperature is at the **Step 1 Start Temp** ± the **Variance**.

 Start program

 Stop program and shut down circulator

 Pause program, press  to restart

Current Cycle: Current Step:  
Total Program Time Remaining (Minutes)


Switching to the Status Display while a ramp program is running shows the status and program time remaining. **Note** If assured soak is enabled the time stops counting down at the end of the step until the desired temperature ± variance is reached. ▲


If a program is not running, highlighting **Ramp** with the circulator running causes the bath temperature to go to the ramp **Start Temp**. The temperature remains there until the ramp is started.

5. To pause the program at the current temperature, highlight the pause program symbol **||** and press .

To restart the program press  again.

6. When the program is complete the circulator either **Shuts Down** or **Maintains** the last setpoint.


When  is pressed after the ramp program is complete the refrigeration, pump and circulator shut off.

**Interfaces** is used to enable/configure the optional serial communications feature.

1. With **Interfaces** highlighted, press  to display the list of parameters.

Highlight the desired parameter and press  to view/change the available options.

**Available options:**

- Serial Type** Off, RS232, RS485 or Analog IO
- Baud** 38400, 19200, 9600, 4800, 2400, 1200, 600 or 300
- Parity** None, Odd or Even
- Data Bits** 8 only
- Stop Bits** 1 or 2
- Address** (Displayed for RS485 only)

Supported protocols: AC, Standard, NC, Namur

See the Appendix for additional information.

**Note** Keypad operation is still available with serial communications enabled. ▲

**Settings - Display Options** is used to view/adjust the circulator's Temperature Units, the Temperature Resolution, the displayed Language, the Display Contrast and the Display Delay.


1. With **Temp. Unit** highlighted press . Use the up/down arrows to highlight the desired temperature scale.

Press .


2. With **Temp. Resolution** highlighted press .

Use the up/down arrows to highlight the desired resolution. Press .


3. With **Language** highlighted press . Use the up/down arrows to highlight the desired language. (Scroll down for additional languages.)

Press .

4. With **Display Contrast** highlighted press .

Press  again and use the up/down arrows to change the contrast. With the desired contrast showing, press  again.


**Note** Holding  for five seconds resets the display contrast to the default level and also brings up the language menu to change, if needed, the displayed language. ▲

5. With **Display Delay** highlighted press  to enable/disable it.

Use the up/down arrows to highlight the time and press  again.

Use the up/down arrows to change the value. Once the desired delay is displayed press .

With **Display Delay** enabled and the Start Display showing, if no arrows are pressed the Start Display changes to the Status Display after the delay expires, see pages 4-2 and 4-3.

**System Messages** is used to view any Warning or Fault messages.

1. With **Messages** highlighted, press  to display the options.

<b>Warnings</b>		^
<b>Faults</b>		
		v
<b>Menu</b>		

**System Run Time** is used to view the circulator (**Unit**) and pump operating hours.

1. With **Run Time** highlighted, press  to display the times.

<b>Unit</b>	<b>xxx hours</b>	^
<b>Pump</b>	<b>xxx hours</b>	
		v
<b>Menu</b>		

**System Configuration** is used to view the circulator's configuration.

1. With **Configuration** highlighted, press  to display the settings.

<b>Head</b>	<b>AC200</b>	
<b>FW</b>	<b>XXXXXXXX.XX</b>	
<b>Checksum</b>	<b>XXXX</b>	
<b>Bath</b>	<b>A10; 115V</b>	
<b>FW</b>	<b>XXXXXXXX.XX</b>	
<b>Menu</b>		

Without a bath connected the display is:

<b>Head</b>	<b>AC200</b>	
<b>FW</b>	<b>XXXXXXXX.XX</b>	
<b>Checksum</b>	<b>XXXX</b>	
<b>Bath</b>	<b>None</b>	
<b>FW</b>	<b>XXXXXXXX.XX</b>	
<b>Menu</b>		

With an invalid bath connected the display is:

<b>Head</b>	<b>AC200</b>	
<b>FW</b>	<b>XXXXXXXX.XX</b>	
<b>Checksum</b>	<b>XXXX</b>	
<b>Bath</b>	<b>Invd xx</b>	
<b>FW</b>	<b>XXXXXXXX.XX</b>	
<b>Menu</b>		

**System - Password/Reset** is used only by a qualified technician. Changing the password enables circulator reset options, the temperature sensor calibration procedure and displays PID values.

1. With **Password/Reset** highlighted, press  to display:

Level	User	^
Password	0	
		v
Menu		

2. Press  and change the number to **1**.

Level	User	^
Password	1	
		v
Menu		

3. Press  to display:

Level	Operator	^
Password	1	
Reset		
Calibration		v
Menu		

4. If desired, highlight **Reset** and press  to display:

Reset user settings	^
Reset PID settings	
Reset both	v
Menu	

**Note** The circulator resets to the **User** mode when it is turned off. The circulator also resets to the **User** mode when the Start/Status Display is displayed continuously for 10 minutes. ▲

Highlight the desired reset option and press .

A confirmation message will appear, press  again.

The circulator will enter the stand by mode.

5. With **Calibration** highlighted, press  to calibrate the internal temperature sensor. The procedure is covered on the next page.

Level	Operator	^
Password	1	
Reset		
Calibration		v
Menu		

6. Scroll down to display **PID Tuning**, see page 4-20.

Password	1	^
Reset		
Calibration		
PID Tuning		v
Menu		

**Note** Ensure the RTA is set to 0 before doing a calibration. ▲

1. To calibrate the temperature sensor highlight **Calibration** and press  to display:

<b>Internal RTD</b>	▲
<b>External RTD</b>	
	▼
<b>Menu</b>	

2. With the desired sensor highlighted, press  to display:

<b>Calibrate</b>	▲
<b>Restore User Cal</b>	
<b>Save User Cal</b>	
<b>Restore Factory Cal</b>	▼
<b>Menu</b>	

3. Press  to display:

<b>Internal RTD</b>	<b>xx.x</b>	▲
<b>High</b>	<b>xx.x</b>	
<b>Low</b>	<b>xx.x</b>	
<b>Calibrate</b>	<b>SP</b>	<b>xx.x</b> ▼
<b>Menu</b>		

4. Highlight the **SP** temperature box and enter either the desired high end or low end setpoint value and press .

<b>Internal RTD</b>	<b>xx.x</b>	▲
<b>High</b>	<b>xx.x</b>	
<b>Low</b>	<b>xx.x</b>	
<b>Calibrate</b>	<b>SP</b>	<b>xx.x</b> ▼
<b>Menu</b>		

5. Depending on which end you are calibrating highlight **High** or **Low**. Once the **Internal RTD** temperature stabilizes for several minutes enter the temperature displayed on your reference thermometer and press .

6. Highlight **Calibrate** and then press  to complete the high end procedure.

<b>Internal RTD</b>	<b>xx.x</b>	▲
<b>High</b>	<b>xx.x</b>	
<b>Low</b>	<b>xx.x</b>	
<b>Calibrate</b>	<b>SP</b>	<b>xx.x</b> ▼
<b>Menu</b>		

<b>Internal RTD</b>	<b>xx.x</b>	▲
<b>High</b>	<b>xx.x</b>	
<b>Low</b>	<b>xx.x</b>	
<b>Calibrate</b>	<b>SP</b>	<b>xx.x</b> ▼
<b>Menu</b>		

7. Repeat for the low end setpoint.

<b>Internal RTD</b>	<b>xx.x</b>	▲
<b>High</b>	<b>xx.x</b>	
<b>Low</b>	<b>xx.x</b>	
<b>Calibrate</b>	<b>SP</b>	<b>xx.x</b> ▼
<b>Menu</b>		

Once complete you can store the calibration into memory by selecting **Save User Cal** and pressing .

You can later restore the same calibration by highlighting **Restore User Cal** and pressing .

Another option is to restore the factory calibration values by highlighting **Restore Factory Cal** and pressing .

1. With **PID Tuning** highlighted, press  to display:

<b>Cool PID</b>	^
<b>Heat PID</b>	
	v
<b>Menu</b>	

2. Highlight the desired PID and press  to display:

<b>P</b>	<b>XX.X</b>	^
<b>I</b>	<b>x.xx</b>	
<b>D</b>	<b>x.xx</b>	
		v
<b>Menu</b>		

3. If required, press  to change the value.

<b>P</b>	<b>XX.X</b>	^
<b>I</b>	<b>x.xx</b>	
<b>D</b>	<b>x.xx</b>	
		v
<b>Menu</b>		


Factory preset values are:

P = 00.6

I = 0.60

D = 0.00

## System - Accessory (AC200 only)


Multifunction Port Pin Out

**Note** Bath circulator accessories are not part of the UL Listings. ▲


**Boost Heater** (pins 5 and 10) is on when enabled and the setpoint is  $>2^{\circ}\text{C}$  above the bath temperature. Boost heater is normally off when the bath temperature is within  $2^{\circ}\text{C}$  of setpoint. **Boost Pump** (pins 4 and 14) is on whenever the option is enabled.

The **Solenoid Valve** and **Auto Refill** are used in conjunction with an *optional* accessory box to select and activate additional features. The box connects to the multifunction port on the rear of the immersion circulator. **Solenoid Valve** for cooling coils (pins 4 and 14) is on when the bath temperature is  $>2^{\circ}\text{C}$  above setpoint or when the heater power drops to 0%. The solenoid is off when the heater power exceeds 80%. **Auto Refill** (pins 5 and 10) is on when the fluid level drops below 47 mm ( $\sim 1\ 7/8$ " ) from the top and is off when the fluid level is at 27 mm ( $\sim 1$ " ) from the top.

1. With **Accessory** highlighted, press  to display:


2. With **Pump/Cool** highlighted press  to display:

Highlight the desired accessory, press  to enable/disable.

3. With **Heater/Auto Refill** highlighted press  to display:


Highlight the desired accessory, press  to enable/disable.

**Refill Timer** is the length of time the auto refill operates if the normal reservoir level is not reached. The range is adjustable from 5 to 600 seconds.

## High Temperature Cutout

HTC (Temperature range varies with type of immersion circulator.)

Ensure the circulator's High Temperature Cutout (HTC) is set to 25°C below the fire point of the selected fluid. ▲

To protect your application, the adjustable High Temperature Cutout (HTC) ensures the heater does not exceed temperatures that can cause serious damage. A temperature sensor is located in the reservoir. A HTC fault occurs when the temperature of the sensor exceeds the set temperature limit.

In the event of a fault the circulator shuts down and displays a fault message, see Section 6. The cause of the fault must be identified and corrected before the circulator can be restarted. A low reservoir fluid level is a primary reason for the HTC to trip.

The HTC is factory preset fully clockwise to the highest possible setting. To set the cutout start the circulator and adjust the setpoint a few degrees higher than the highest desired fluid temperature. Allow the circulator to stabilize at the temperature setpoint. Using a small diameter (less than 1/8") plastic/metal rod, turn the arrow counterclockwise until the circulator shuts down and the fault message appears. Press  to clear the message.

**Note** If using plastic rod ensure the fluid temperature is low enough to prevent possibly melting the plastic. ▲

Before you can restart the circulator it has to cool down a few degrees. To restart the circulator press the black reset ring surrounding the red dial - and then press  again. If Auto Restart is enabled the circulator restarts, if disabled use the Start Up procedure.

**Note:** We recommend periodically rechecking operation or if the circulator is moved. ▲

## Decommissioning/ Disposal

Decommissioning prepares equipment for safe and secure transportation.

**Laboratory Grade Ethylene glycol (EG) is poisonous and flammable. Before disposing refer to the manufacturer's most current SDS for handling precautions. ▲**

**Decommissioning must be performed only by qualified dealer using certified equipment. All prevailing regulations must be followed. ▲**

A45HC units utilize flammable refrigerant (R290), risk of fire or explosion. Dispose of properly in accordance with federal or local regulations.

Consider decommissioning the circulator when:

- It fails to maintain desired specifications
- It no longer meets safety standards
- It is beyond repair for its age and worth

Refrigerant (R134A) and oil (POE or Ester) used in all other refrigerated baths must be recovered from equipment before disposal.

**Note** Keep in mind any impact your application may have had on the circulator. ▲

Direct questions about circulating bath decommissioning or disposal to our Sales, Service and Customer Support.

**Handling and disposal should be done in accordance with the manufacturer's specification and/or the SDS for the material used. ▲**

## Storage

**If the circulator is to be transported and/or stored in cold temperatures it needs to be drained and then flushed with a 50/50 laboratory grade glycol/water mixture. ▲**

The circulator can be stored for up to 90 days inside the temperature range of -25°C to +60°C (-13°F to +140°F).

If necessary, when removing the circulator from storage allow it time to warm up and dry out in order to prevent any condensation issues.

## Section 5 Preventive Maintenance

Disconnect the power cord prior to performing any maintenance. ▲

Handle the circulator with care. Sudden jolts or drops can damage its components. ▲

### Cleaning

After time, the circulator's stainless steel surfaces may show spots and become tarnished. Normal stainless steel cleaners can be used.

Clean the bath vessel and built-in components at least every time the reservoir fluid is changed. Use water and a soft cloth.

Do not use scouring powder. ▲

The inside of the bath must be kept clean in order to ensure a long service life. Quickly remove substances containing acidic or alkaline substances and metal shavings as they could harm the surfaces causing corrosion. If corrosion (e.g., small rust marks) occurs in spite of this, cleaning with stainless steel caustic agents has proved to be suitable. Apply these substances according to the manufacturer's recommendations.

### Condenser Fins

Do not use any substances which contain solvents. ▲

In order to maintain the cooling capacity of a refrigerated bath, clean the fins two to four times per year, depending on the operating environment.

Switch off the circulator and unplug the power cord.

For ARCTIC A45HC

1. Remove the condenser panel.
2. Clean fins with brush or similar tool.
3. Replace the panel.

For all other refrigerated baths:

Clean the fins with compressed air.

For extreme soiling a qualified technician will need to remove the cooling compressor casing.

### Tubing


Inspect and tighten the tubing and clamps daily.

## Grounding Strap and Nut


Sahara baths, and stand-alone circulators, require immersion circulator grounding.

If the circulator is removed for any reason you must ensure the nut and washer secures to the grounding strap located on the rear of the tank top when the circulator is reinstalled. ▲

## Testing the Safety Features

### High temperature protection

Using a small diameter (less than 1/8") plastic/metal rod, turn the arrow to the desired temperature.

**Note** If using plastic rod ensure the fluid temperature is low enough to prevent possibly melting the plastic. ▲

Set a cut-off temperature that is lower than the desired setpoint temperature.

Switch on the circulator and ensure it shuts down at the set cut-off temperature.

After the circulator cools down clear the HTC error message by pressing the black ring surrounding the red dial - and then press . If Auto Restart is enabled the circulator restarts, if disabled use the Start Up procedure.

If the circulator did not shut down have it checked by a qualified technician.

Reset the safety to the desired temperature.

### Low liquid level protection

With the circulator on, use a screwdriver and slowly push down on each level sensor, one at a time, until an error message appears. See Section 6 for details on error messages.

If not, have the circulator checked by a qualified technician.

## Error Displays

### Section 6 Troubleshooting

The circulator can display error messages and, if enabled, sound an alarm. Error messages are cleared by pressing enter (↵). Restart the circulator once the cause of the error message is identified and corrected. If the cause was not corrected the error code will reappear, contact our Sales, Service and Customer Support. If **Auto start** is enabled the circulator will restart, if disabled use the Start Up procedure.


**FAULT:**  
**High Temperature**  
**PRESS ENTER**  
**to clear message**

Error Display (Typical)

#### Fault Displays

The heating element, pump and, if applicable, refrigeration shut down with a fault. A fault also sounds the alarm, if enabled.

Message	Cause	Actions
<b>Bath not found</b>	<ul style="list-style-type: none"> <li>a loss of communication between refrigerated bath and circulator</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (↵)</li> <li>check cable connections on the rear of the bath and circulator</li> </ul>
<b>High Fixed Temp.</b>	<ul style="list-style-type: none"> <li>circulator's nonadjustable high temperature protection limit exceeded</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (↵)</li> <li>check fluid selection</li> <li>check environmental conditions</li> </ul>
<b>High Temperature Refrigeration</b>	<ul style="list-style-type: none"> <li>high refrigeration temperature</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (↵) check voltage supply</li> <li>the refrigeration may need servicing, contact us</li> </ul>
<b>HPC High Press. Cutout</b>	<ul style="list-style-type: none"> <li>the high refrigeration pressure cutout activated</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (↵)</li> <li>check for obstructions to air flow</li> <li>the refrigeration may need servicing, contact us</li> </ul>

Message	Cause	Actions
<b>HTC-LLC</b>	<ul style="list-style-type: none"> <li>if the fluid level is greater than approximately 57 mm (2 1/4") below the reservoir top it is a low level fault.</li> <li>if not, high temperature protection limit exceeded</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (↵)</li> <li>if required, check for leaks and fill fluid to proper level</li> <li>on smaller baths, 2 - 4 liters, if the fluid is near the minimum level reducing the pump speed may eliminate nuisance messages.</li> <li>if fluid is not the issue, allow circulator to cool down to at least 10°C below the HTC setting</li> <li>turn the red knob on the HTC fully clockwise</li> <li>press the HTC's black reset ring</li> <li>press enter (↵) again</li> <li>reset HTC to desired setting, see Section 4</li> <li>if message reappears recycle power to circulator and repeat the procedure</li> <li>increasing the pump speed to maximum when using fluids with viscosity greater than 50Cst may eliminate nuisance messages</li> </ul> 
<b>High Level</b>	<ul style="list-style-type: none"> <li>reservoir fluid level too high for safe operation</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (↵)</li> <li>check fluid level, drain excess fluid if required</li> <li>verify optional auto refill operation</li> </ul>
<b>LLC Low Level Cutout</b>	<ul style="list-style-type: none"> <li>reservoir fluid level too low for safe operation</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (↵)</li> <li>check fluid level</li> <li>check for leaks</li> </ul>
<b>Low Fixed Temp.</b>	<ul style="list-style-type: none"> <li>circulator's nonadjustable low temperature protection limit exceeded</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (↵)</li> <li>check fluid selection</li> </ul>

Message	Cause	Actions
<b>Low Temperature</b>	<ul style="list-style-type: none"> <li>adjustable high temperature protection limit exceeded</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>check limit setting</li> <li>check fluid selection</li> </ul>
<b>Motor Fault</b>	<ul style="list-style-type: none"> <li>pump is not turning due to possible motor overload</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>it can take over 10 minutes for the motor temperature to get low enough before the circulator can be restarted</li> <li>if the motor does not start, cycle the power</li> <li>if multiple attempts do not clear the issue the motor probably needs to be replaced, contact us</li> </ul>
<b>MOL Motor Overload</b>	<ul style="list-style-type: none"> <li>high motor overload temperature</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>allow circulator to cool down</li> </ul>
<b>Open RTD1 Internal</b>	<ul style="list-style-type: none"> <li>open internal temperature sensor</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>contact us</li> </ul>
<b>Open RTD2 External</b>	<ul style="list-style-type: none"> <li>open external temperature sensor</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>contact us</li> </ul>
<b>Refill Error</b>	<ul style="list-style-type: none"> <li>Reservoir did not fill within the user specified auto refill time</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>check water source is turned on</li> <li>check multifunction port connection</li> <li>check time specified</li> <li>check for leaks</li> </ul>
<b>Shorted RTD1 Internal</b>	<ul style="list-style-type: none"> <li>shorted internal temperature sensor</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>contact us</li> </ul>
<b>Shorted RTD2 External</b>	<ul style="list-style-type: none"> <li>shorted external temperature sensor</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>contact us</li> </ul>

### Warning Displays

The circulator will continue to run with a warning. A warning also sounds the alarm, if enabled.

Message	Cause	Actions
<b>Bad Calibration</b>	<ul style="list-style-type: none"> <li>bad temperature probe calibration</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>redo calibration</li> </ul>
<b>High Temperature</b>	<ul style="list-style-type: none"> <li>adjustable high temperature protection limit exceeded</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>check limit setting</li> <li>check fluid selection</li> </ul>
<b>Low Level</b>	<ul style="list-style-type: none"> <li>reservoir fluid level too low for safe operation</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>check fluid level</li> <li>on smaller baths, 2 - 4 liters, if the fluid is near the minimum level reducing the pump speed may eliminate nuisance messages.</li> </ul>
<b>Low Temperature</b>	<ul style="list-style-type: none"> <li>adjustable low temperature protection limit exceeded</li> </ul>	<ul style="list-style-type: none"> <li>clear message by pressing enter (  )</li> <li>check limit setting</li> <li>check fluid selection</li> </ul>

## Checklist

### Circulator will not start

Check the display for error codes, see Error Codes in this section.

Ensure the circuit protector(s) is in the on ( I ) position.

Make sure supply voltage is connected and matches the circulator's nameplate rating  $\pm 10\%$ .

### No display

Recycle the circuit protector on the rear of the circulator.

### Display reads 239°C

External temperature probe selected but no probe attached to circulator.

Attach external probe or select internal temperature sensor.

### Circulator will not circulate process fluid

Check the reservoir level. Fill, if necessary.

Check the application for restrictions in the cooling lines.

The pump motor overloaded. The pump's internal overtemperature overcurrent device will shut off the pump causing the flow to stop. This can be caused by low fluid, debris in system, operating circulator in a high ambient temperature condition or excessively confined space. Allow time for the motor to cool down.

Ensure supply voltage matches the circulator's nameplate rating  $\pm 10\%$ .

### Inadequate temperature control

Verify the setpoint.

For refrigerated baths, ensure the condenser is free of dust and debris.

Check the fluid concentration.

Ensure circulator installation complies with the site requirements in Section 3.

Ensure supply voltage matches the nameplate rating  $\pm 10\%$ .

If the temperature continues to rise, ensure your application's heat load does not exceed the rated specifications.

Enter the controller menu and ensure the ENERGY SAVER mode is on in order for the system to maintain a stable temperature.

Check for high thermal gradients (e.g., the application load is being turned on and off or rapidly changing).

### **Inadequate pump pressure**

Ensure any user installed in-line valves are in the desired position.

Ensure the circulator's process fluid outlet is connected to the application's fluid inlet and not the outlet, see Section 3.

Ensure all connections are secure and that the proper sealant/lubricant for the fitting material is used.

Keep the distance between the circulator and the application as short as possible.

Ensure tubing is straight and without bends. If diameter reductions are required, make them at the inlet and outlet of your application, not at the circulator.

### **Circulator shuts down**

Ensure  wasn't accidentally pressed.

Ensure the circuit protector(s) is in the on (I) position.

Check the display for error codes.

Make sure supply voltage is connected and matches the nameplate rating  $\pm 10\%$ .

Restart the circulator.

### **Circulator fan not running**

The fan comes on when internal head temp is greater than 40°C and shuts off when the temperature drops below 35°C.

### **USB Driver Not Recognized**

If your operating system does not automatically recognize the optional driver log on to:

<http://www.ftdichip.com/FTDrivers.htm>

for instructions.

Please contact Thermo Fisher Scientific Sales Service and Customer Support if you need any additional information, see inside cover for contact instructions.

# Appendix AC Serial Communications Protocol

Serial communication is accomplished either through the optional 9-pin Serial Communications Box or through the USB port on the immersion circulator. If your operating system does not automatically recognize the optional driver log on to: <http://www.ftdichip.com/FTDrivers.htm> for instructions.

**Note** This appendix assumes you have a basic understanding of communications protocols. Information on the NC, STANDARD and NAMUR protocols is available upon request.

**Note** Keypad operation is still available with serial communications enabled.

**Note** NC protocol is required to use RS485 device addressing.

All commands must be entered in the exact format shown in the tables on the following pages. The tables show all commands available, their format and responses. Controller responses are either the requested data or an error message. The controller response must be received before the host sends the next command.

The host sends a command embedded in a single communications packet, then waits for the controller's response. If the command is not understood, the controller responds with an error command. Otherwise, the controller responds with the requested data.

Commands are not case sensitive. Upper or lower case letters may be used. Commands are listed in the Commands Table, error responses are given in the Errors Table, and symbols are shown in the Key Table.

Key	
Symbol	Meaning
[B]	A binary value 0 or 1 (0 = Off, FALSE or Disable(d); 1 = On, TRUE or Enable(d)).
[CR]	Carriage return – used as the termination character.
[U]	Text representing the units associated with a value.
[V]	A value that can be requested in a read command or sent as part of a set command.
[V <sub>MAX</sub> ]	Maximum allowed value. Part of error message when set value is too high.
[V <sub>MIN</sub> ]	Minimum allowed value. Part of error message when set value is too low.

Value: Read commands return analog [V] or bit [B] values or settings, while set commands send analog or bit settings. Read commands return values in the same displayed precision Set command messages missing the space character between the command and the setting will be rejected, as the user's intent is unclear.

Units: A read command returning an analog [V] value or setting, will include the units [U] associated with that value or setting. A set command sending an analog value will not include the units. The units returned by the complementary read command are assumed.

Termination character: A carriage return [CR] is used to terminate command and response messages. (Typically the "Enter" key on the keyboard.)

**Note** The inter-character timeout (time between transmitted characters) is set to 30 seconds. Exceeding the timeout will clear the receiver buffer and require the message to be retransmitted.

**Note** Special characters (backspace, delete, insert, etc.) are not recognized and generate error responses.

## Commands Table:

<b>Commands</b>		<i>All messages from master and slave are terminated with a carriage return [CR]</i>	
Command Description	<i>Notes</i>	Master Sends	Sample Slave Response <i>(echo off)</i> <span style="float: right;"><i>Alternate units</i></span>
Read Temperature	<i>Internal</i>	RT	[V]C <span style="float: right;">F K</span>
Read Temperature 2	<i>External</i>	RT2	[V]C <span style="float: right;">F K</span>
Read Displayed Setpoint		RS	[V]C <span style="float: right;">F K</span>
Read Setpoint X (X = 1 to 5)		RSX	[V]C <span style="float: right;">F K</span>
Read High Temperature Fault		RHTF	[V]C <span style="float: right;">F K</span>
Read High Temperature Warn		RHTW	[V]C <span style="float: right;">F K</span>
Read Low Temperature Fault		RLTF	[V]C <span style="float: right;">F K</span>
Read Low Temperature Warn		RLTW	[V]C <span style="float: right;">F K</span>
Read Proportional Heat Band Setting		RPH	[V]%
Read Proportional Cool Band Setting		RPC	[V]%
Read Integral Heat Band Setting		RIH	[V]Repeats per minute
Read Integral Cool Band Setting		RIC	[V]Repeats per minute
Read Derivative Heat Band Setting		RDH	[V]Minutes
Read Derivative Cool Band Setting		RDC	[V]Minutes
Read Temperature Precision		RTP	[V]
Read Temperature Units		RTU	[V} <span style="float: right;">C,F,K</span>
Read Unit On		RO	[B]
Read External Probe Enabled		RE	[B]
Read Auto Restart Enabled		RAR	[B]
Read Energy Saving Mode		REN	[B]
Read Time		RCK	hh:mm:ss
Read Date		RDT	mm/dd/yyyy or dd/mm/yyyy
Read Date Format		RDF	MM/DD/YYYY or
Read Ramp Status		RRS	<i>Stopped, Running, Paused</i>

<b>Commands</b>		<i>All messages from master and slave are terminated with a carriage return [CR]</i>	
Command Description	Notes	Master Sends	Sample Slave Response <i>(echo off)</i> <span style="float: right;"><i>Alternate units</i></span>
Read Firmware Version		RVER	[V]
Read Firmware Checksum		RSUM	[V]
Set Displayed Setpoint		SS [V]	OK
Set Setpoint X (X = 1 to 5)	See Note 1	SSX [V]	OK
Select Setpoint X (X = 1 to 5)	See Note 1	SSX	OK
Set High Temperature Fault		SHTF [V]	OK
Set High Temperature Warning		SHTW [V]	OK
Set Low Temperature Fault		SLTF [V]	OK
Set Low Temperature Warning		SLTW [V]	OK
Set Proportional Heat Band Setting		SPH [V]	OK
Set Proportional Cool Band Setting		SPC [V]	OK
Set Integral Heat Band Setting		SIH [V]	OK
Set Integral Cool Band Setting		SIC [V]	OK
Set Derivative Heat Band Setting		SDH [V]	OK
Set Derivative Cool Band Setting		SDC [V]	OK
Set Temperature Resolution		STR [V]	OK
Set Temperature Units		STU [V] C,F,K	OK
Set Unit On Status		SO [B]	OK
Set External Probe On Status		SE [B]	OK
Set Auto Restart Enabled		SAR [B]	OK
Set Energy Saving Mode		SEN [V]	OK
Set Pump Speed		SPS [V] L,M,H	OK
Set Ramp Number		SRN [V]	OK

Note 1: SSX is a content dependent command. If SSX is immediately followed by a space and a value, the value is assigned to setpoint x. For example, if setpoint 2 is currently 37° and SS2 45<cr> is issued, setpoint 2 is changed to 45°. If SSX is immediately followed by a carriage return the controlling setpoint is changed to X. For example, if the controlling setpoint is 3 and SS<cr> is issued, the controlling setpoint is changed to setpoint 2.

<b>Commands</b>		<i>All messages from master and slave are terminated with a carriage return [CR]</i>	
Command Description	Notes	Master Sends	Sample Slave Response <i>(echo off) Alternate units</i>
Set Ramp Program	<i>this will load defaults for all steps declared by V1</i>	SRP [v1,v2,v3,v4]  V1 #of steps V2 variance V3 cycles V4 start temp	OK  Note: If unit allows more than one ramp the ramp number must first be set.
Set Ramp Step		SRS[v1,v2,v3]  V1 step # V2 end temp V3 duration	OK  Note: If unit allows more than one ramp the ramp number must first be set and SRP defined.
Set Ramp On Status		SRO [V] S,E,P  Start, End, Pause/resume	OK  Note: (P)pause toggles Pause/resume

Errors Table:

<b>Errors</b>	
Error Description	Notes Slave Responds
<i>Not defined, not implemented or incorrectly formatted</i>	? Unsupported command
<i>Extra characters...</i>	? Format error
<i>Set value too high</i>	? Maximum allowed is $[V_{MAX}]$
<i>Set value too low</i>	? Minimum allowed is $[V_{MIN}]$
<i>Argument to binary set command not 0 or 1</i>	? Value must be 0 or 1
<i>Set command attempted while in read only mode</i>	? Mode is read only
<i>Set command failed (e.g. SO 1 with low level)</i>	? Failed

**RUFS** Read Unit Fault Status

This command returns 5 values. These are decimal representations of hexadecimal values. Each individual bit of the value represents a different warning, fault or status.

decimal	hex	B7	B6	B5	B4	B3	B2	B1	B0
1	1	0	0	0	0	0	0	0	1
2	2	0	0	0	0	0	0	1	0
4	4	0	0	0	0	0	1	0	0
8	8	0	0	0	0	1	0	0	0
16	10	0	0	0	1	0	0	0	0
32	20	0	0	1	0	0	0	0	0
64	40	0	1	0	0	0	0	0	0
128	80	1	0	0	0	0	0	0	0

Value	Description of bits
$V_1$	B0 - b5 unused B6 rtd1 shorted B7 rtd1 open
$V_2$	B0 HTC fault B1 high RA temperature fault B2 - b7 unused
$V_3$	B0 low level warn B1 lo temperature warn B2 high temperature warn B3 low level fault B4 lo temperature fault B5 high temperature fault B6 low temperature fixed fault B7 high temperature fixed fault
$V_4$	B0 PWM heat duty cycle > 0 B1 compressor On/Off B2 Pump On status B3 Unit On status B4 Unit Stopping B5 Unit fault status B6 unused B7 Beeper On status
$V_5$	B0 Pump speed fault B1 MOL fault B2 HPC fault B3 Cool Icon On steady (unit is cooling at max capacity) B4 Cool Icon flashing (unit is B5 Heat Icon On steady B6 Heat Icon flashing B7 External sensor controlling

Refer to Key table on page 1 for explanation of symbols and their meanings.

Examples:

Read Temperature:

Host

R	T		<b>CR</b>
Command			[CR]

Controller:

2	0	.	0	C	<b>CR</b>
[V]			[U]	[CR]	

Set Setpoint:

Host

S	S		2	0	<b>CR</b>
Command			[V]	[CR]	

Controller:

O	K	<b>CR</b>
Command Accepted		[CR]

Read Temperature 2:

Host:

R	T	2	<b>CR</b>		
2	0	.	0	C	[CR]

Controller:

Set Setpoint to -22°C when minimum allowed is -20°C: Minimum allowed is  $[V_{MIN}]$

Host:

S	S		-	2	2	<b>CR</b>															
?		M	i	n	i	m	u	m	a	l	l	o	w	e	d	i	s	-	2	0	<b>CR</b>

Controller:

# Warranty

Thermo Fisher Scientific warrants for 36 months from date of shipment the Thermo Scientific ADVANCED series of Immersion Circulators, ARCTIC refrigerated bath circulators, and SAHARA heated bath circulators according to the following terms.

Any part of the unit manufactured or supplied by Thermo Fisher Scientific and found in the reasonable judgment of Thermo Fisher to be defective in material or workmanship will be repaired at an authorized Thermo Fisher Repair Depot without charge for parts or labor. The unit, including any defective part must be returned to an authorized Thermo Fisher Repair Depot within the warranty period. The expense of returning the unit to the authorized Thermo Fisher Repair Depot for warranty service will be paid for by the buyer. Our responsibility in respect to warranty claims is limited to performing the required repairs or replacements, and no claim of breach of warranty shall be cause for cancellation or rescission of the contract of sales of any unit. With respect to units that qualify for field service repairs, Thermo Fisher Scientific's responsibility is limited to the component parts necessary for the repair and the labor that is required on site to perform the repair. Any travel labor or mileage charges are the financial responsibility of the buyer.

The buyer shall be responsible for any evaluation or warranty service call (including labor charges) if no defects are found with the Thermo Scientific product.

This warranty does not cover any unit that has been subject to misuse, neglect, or accident. This warranty does not apply to any damage to the unit that is the result of improper installation or maintenance, or to any unit that has been operated or maintained in any way contrary to the operating or maintenance instructions specified in this Instruction and Operation Manual. This warranty does not cover any unit that has been altered or modified so as to change its intended use.

In addition, this warranty does not extend to repairs made by the use of parts, accessories, or fluids which are either incompatible with the unit or adversely affect its operation, performance, or durability.

Thermo Fisher Scientific reserves the right to change or improve the design of any unit without assuming any obligation to modify any unit previously manufactured.

THE FOREGOING EXPRESS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO WARRANTIES OR MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

OUR OBLIGATION UNDER THIS WARRANTY IS STRICTLY AND EXCLUSIVELY LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE COMPONENT PARTS AND Thermo Fisher Scientific DOES NOT ASSUME OR AUTHORIZE ANYONE TO ASSUME FOR IT ANY OTHER OBLIGATION.

Thermo Fisher Scientific ASSUMES NO RESPONSIBILITY FOR INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES INCLUDING, BUT NOT LIMITED TO LOSS OR DAMAGE TO PROPERTY, LOSS OF PROFITS OR REVENUE, LOSS OF THE UNIT, LOSS OF TIME, OR INCONVENIENCE.

This warranty applies to units sold in the United States. Any units sold elsewhere are warranted by the affiliated marketing company of Thermo Fisher Scientific. This warranty and all matters arising pursuant to it shall be governed by the law of the State of New Hampshire, United States. All legal actions brought in relation hereto shall be filed in the appropriate state or federal courts in New Hampshire, unless waived by Thermo Fisher Scientific.

Thermo Fisher Scientific  
168 Third Avenue  
Waltham, Massachusetts 02451  
United States

[www.thermofisher.com](http://www.thermofisher.com)