

ermoscientifi

# TSX Series high-performance undercounter refrigerators

The only undercounter lab-grade refrigerator with redundant cooling

- Low energy usage
- Low noise
- High storage capacity
- Universal voltage





Thermo Fisher SCIENTIFIC

# Introducing the Thermo Scientific™ TSX505

### Powered by V-Drive featuring STeM Technology

- Redundant cooling technology protects your storage while delivering 37% less energy, 20X less noise and twice the storage capacity\*
- Variable speed control using thermoelectric devices in place of compressors, provides ideal conditions for the most demanding applications
- Universal voltage (100-240v, 50/60hz) allows placement throughout the facility and globally just by connecting the region specific power cord.



## Designed for the following industries



\*Internal data based on comparisons to Thermo Scientific model

# How the TSX505 works

### Compressor-free V-drive featuring STeM (synchronized temperature management)

- Dual refrigeration system combines the benefits of cold wall and fan-assisted recovery
- · Actively monitors the cabinet to maintain the desired temperature conditions
- Constant monitoring and adjustment during periods of door openings keeps temperatures stable without temperature fluctuation found in models with compressors that cycle on/off
- This tight internal air circulation eliminates the risk of freezing which can ruin materials like vaccines



The TSX505's refrigerant loop absorbs and transfers heat energy from the outside into a heat pump in the back of the refrigerator. There, heat rejection coils transfer heat energy from the heat pump back into the ambient environment.

### System Comparison Chart

In the chart, the gray line represents the additive effect of using cold wall (blue line) and forced air (orange line), with the benefit that if one of those systems should experience an issue, you know the other system can keep the unit going until service arrives.



# Designed to save energy, reduce noise and meet sustainability protocols

#### As shown, the TSX505 uses up to 37% less energy with an annual cost reduction of 39%\*

Energy Usage				
Refrigerator Model	Annual Power Usage (kW)	Annual Run Time (hours)	Annual Energy Usage (kWh)	Additional Energy Used
TSX505	0.11	24	2.7	0%
REL404	0.18	24	4.30	37.2%

Annual Cost to Run								
Refrigerator Mode	Cycle/Day	Usage Hours/Day	Usage Hours/Year	Energy Cost (\$/kWh)**	Annual Cost	Additional Cost		
TSX505	1	24	8,760	\$0.1019	\$98.19	0%		
REL404	1	24	8,760	\$0.1019	\$160.68	30%		

#### The sound of quiet

Whisper-quiet operation is the unseen feature of the TSX505. As shown, the 35 dBA sound output is less than an average library, enabling the refrigerator to reside comfortably in work station areas.

#### The SNAP/F-Gas compliant choice

The TSX505 is designed to meet the standards established by the U.S. Environmental Protection Agency's Significant New Alternatives Policy (SNAP) program and the European Union's F-Gas Regulation (517/2014). SNAP evaluates refrigerants that are potentially harmful to the environment. The TSX505 has been designed with refrigerants that are in compliance with SNAP.

\*Internal data based on comparisons to Thermo Scientific model REL404A. \*\*https://www.eia.gov/electricity/monthly/epm\_table\_grapher.cfm?t=epmt\_5\_6\_a. Use data from the most current month reported and use the energy rate for the "US Total" in the "Commercial" sector.

## Decibel scale of common sounds

Noise source	Decibel level (dBA)
Freeway traffic, vacuum cleaner	70
Conversation in restaurant, office, background music	60
Library	40
TSX Series undercounter refrigerat	or 35†
Quiet rural area	30
Ductling looved	20
Rustling leaves	20

† Steady state 35 dBa; start-up 42 dBa

# Designed to hold more within a smaller footprint



As shown in these two side views, the TSX505 does not include internal protrusions such as overhead fans. This allows you to store more inside its interior. Also, with the V-drive and STeM cooling, temperatures are consistent at all shelf levels, making the TSX505 a smart choice for vaccines, reagents, media and other sensitive material.

# 495 BOXES

TSX Series undercounter

- 45 boxes per level (9x5 pattern)
- 225 boxes per shelf (5 levels per shelf, 2 shelves per unit)
- +45 boxes (1 additional level above top shelf)

#### Total number of vaccine boxes: 495 boxes (225 + 225 + 45)\*\*

\*Thermo Scientific Revco, model REL404 \*\*MMR vaccine box: 3.5"W x 1.75"H x 1.625"D

## Alternative model\*



- 36 boxes per level (9x4 pattern)
- 108 boxes per shelf (3 levels per shelf, 2 shelves per unit)

Total number of vaccine boxes: 216 boxes (108 + 108)

#### Easy to operate user interface

- Password protected
- Icon based
- Adjustable lighting
- Built-in datalogging
- USB port
- Ergonomic tilted design to help prevent neck strain
- Electronic lock is standard



## thermo scientific

#### Specifications and ordering information

TSX Series high-performance undercounter refrigerator													
Model No.	Doors	Settable Range	Energy Consumption	Capacity cu. ft. (liters)	Electrical	Plug Type	Shelves	Defrost	Sound Level	Certifications	Interior Dimensions: (D x W x H)	Exterior Dimensions: (D x W x H)	Shipping Weight Ibs. (kg)
TSX505SA	1, Solid	+3°C to	1.7 kWh/day	156 Liters - / 5.5	100 - 240v, 50/60hz*	NEMA 5-15	3	Auto	35 dBa 3	UL, CE, REACH, RoHS	20.5 x 19.5 x 23 in. (520 x 495 x 584 mm)	26 x 23.6 x 31.8 in. (660.4 x 599.4 x 807.7 mm)	155 (71)
TSX505GA	1, Glass	+7°C	2.4 kWh/day	Cubic Feet									147 (66.7)

\*For alternate power cords, review accessory list below.

#### **Customer installed accessories**

For use with refrigerator size	5.5 cu. ft. (156L)			
Description	Cat. No.			
Rolling Wire Basket	WB05			
Additional wire shelf	WS05			
Set of casters with mounting hardware	CS05			
Power Cord 115V	100056H01			
Power Cord, D NEMA 6-15P 230V	100091H01			
Power Cord, V EU Schuko 230V	100089H01			
Power Cord, C China 220V/50Hz	100092H01			
Power Cord, W UK Type G-10A 23	100090H01			

#### Thermo Scientific<sup>™</sup> Snappies<sup>™</sup> Breast Milk Containers



Thermo Scientific<sup>™</sup> Nalgene<sup>™</sup> Bottles and Containers



Media and Reagents

**Thermo Fisher** S C I E N T I F I C

Visit **thermofisher.com** for accessories, wireless monitoring, consumables, media, reagents and more.

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for the customers specific use or application. © 2017 - 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. The SNAP program is run by the U.S. Environmental Protection Agency. The Energy Star trademark belongs to the U.S. Environmental Protection Agency and the U.S. Department of Energy. BRCSTSXUC 0620