




MSAA111A1

\$19

Product Images



Description

Diffuser for “AM” Dino-Lite [Basic & Premier Series](#), excluding “Z” Polarized models. Proper installation of the diffuser: Place the diffuser into the clear plastic nozzle. Apply slight pressure until diffuser is secure. Please be caution of misaligned Led lights. Proper removal of the diffuser: Grip the clear plastic nozzle of the microscope with one hand and with the other hand gently twists and removes the diffuser.

Additional Information

Warranty Period	1 year
Manufactured in	China
Service & Support	U.S. (www.dinolite.us) English & Spanish
Regulatory Approval	ROHS
Dimensions	29mm(D) x 23mm(L)
Weight	3g
Material	Plastic
Compatible with	“AM” Basic & Premier Series, Excluding “Z” Models
Package Includes	Diffuser, Instruction Guide
Package Dimensions	7.6cm(L) x 5.1cm(W) x 3.2(H)
Package Weight	20g

